



TECHNICAL DATA SHEET

Sterile Fabric Adhesive Bandage

Items # AF0310 to AFK0315-B

Intended Use

Sterile fabric adhesive bandage is used to cover and protect wounds or to hold together the skin edges of a wound.

Characteristic	Specification				
Description	Die-cut fabric adhesive bandage with absorbent non-stick pad in the center and adhesive side completely covered with two non-stick paper strips.				
Material	Fabric: 100% Polyester				
Color	Skin				
Dimension	Item Code	Width (cm)		Length (cm)	
	AF0310, AF0310-B, AF0310-H	2,5		7,6	
	AF0328	7,3			
	AF0334	1,9			
	AF0402	5,1		10,2	
	AFF0217, AFF0317-B	4,2		5,1	
	AFK0315, AFK0315-B			7,6	
Packaging	Item Code	Bandage/ Pouch	Pouches/ Box	Boxes/ Case	Bandages/ Case
	AF0310, AF0310-H, AFK0315, AF0334, AFF0217	1	100	12	1200
	AF0402		50	24	
	AF0328			12	600
	Item Code (Bulk)	Bandage/ Pouch	Pouches/ Case		Bandages/ Case
	AFK0315-B	1	1800		1800
	AFF0317-B		2400		2400
	Item Code (Bulk)	Bandages/Tray	Trays/ Case		Bandages/ Case
	AF0310-B	130	10		1300

Item Code	Description	Color	Specific Attributes
AF0310	AMD-Ritmed Fabric Adhesive Bandages – 1" x 3" (2,5 cm x 7,6 cm), 1 Bandage/Pouch, 100 Pouches/Box, 12 Boxes/Case, 1200 Bandages/Case	Skin	Sterile
AF0310-B	AMD-Ritmed Fabric Adhesive Bandages – 1" x 3" (2,5 cm x 7,6 cm), 130 Bandages/Tray, 10 Trays/Case, 1300 Bandages/Case	Skin	Sterile, Bulk
AF0310-H	AMD-Ritmed Heavy Weight Fabric Bandages – 1" x 3" (2,5 cm x 7,6 cm), 1 Bandage/Pouch, 100 Pouches/Box, 12 Boxes/Case, 1200 Bandages/Case	Skin	Sterile
AF0328	AMD-Ritmed Flexible Fabric 4 Wing Adhesive Bandage – 2 7/8" x 3" (7,3 cm x 7,6 cm), 1 Bandage/Pouch, 50 Pouches/Box, 12 Boxes/Case, 600 Bandages/Case	Skin	Sterile
AF0334	AMD-Ritmed Fabric Adhesive Bandages – 3/4" x 3" (1,9 cm x 7,6 cm), 1 Bandage/Pouch, 100 Pouches/Box, 12 Boxes/Case, 1200 Bandages/Case	Skin	Sterile
AF0402	AMD-Ritmed Fabric Jumbo Adhesive Bandage – 2" x 4" (5,1 cm x 10,2 cm), 1 Bandage/Pouch, 50 Pouches/Box, 24 Boxes/Case, 1200 Bandages/Case	Skin	Sterile



TECHNICAL DATA SHEET
Sterile Fabric Adhesive Bandage
Items # AF0310 to AFK0315-B

Item Code	Description	Color	Specific Attributes
AFF0217	AMD-Ritmed Fabric Fingertip Bandages – 1 3/4" x 2" (4,4 cm x 5,1 cm), 1 Bandage/Pouch, 100 Pouches/Box, 12 Boxes/Case, 1200 Bandages/Case	Skin	Sterile
AFF0317-B	AMD-Ritmed Fabric Fingertip Adhesive Bandages – 1 3/4" x 2" (4,4 cm x 5,1 cm), 1 Bandage/Pouch, 2400 Pouches/Case, 2400 Bandages /Case	Skin	Sterile, Bulk
AFK0315	AMD-Ritmed Fabric Knuckle Adhesive Bandages – 1 1/2" x 3" (3,8 cm x 7,6 cm), 1 Bandage/Pouch, 100 Pouches/Box, 12 Boxes/Case, 1200 Bandages/Case	Skin	Sterile
AFK0315-B	AMD-Ritmed Fabric Knuckle Adhesive Bandages – 1 1/2" x 3" (3,8 cm x 7,6 cm), 1 Bandage/Pouch, 1800 Pouches/Case, 1800 Bandages /Case	Skin	Sterile Bulk

Regulatory

FDA Class 1, FDA Product Code KGX

Health Canada Class 1

Storage Conditions

Store under ambient temperature and dry conditions.