

TECHNICAL DATA SHEET

Trauma Gauze & Nasal Dressing Items # AHM0022 to AHMDG20

Intended Use

Used in Health Care facilities for management of nasal orifice and topical wounds and temporary control of external surface bleeding such as minor cuts and minor lacerations.

Characteristic	Specification				
Description	Item code AHM0216 (Nasal Dressing): Hard core covered with a gauze sleeve. The sleeve is closed at one end by a double string stitched through the hard cord and at the other end by being folded inward and fixed. Other item codes (Trauma Gauze) and item code AHM0126 gauze sleeve: Soft, glossy, slightly stretchable woven material with cut and/or weaved edges. Item code AHMDG20 is folded in 4 and opens up as a strip.				
Material Composition	Gauze	Carboxymethyl cellulose			
	Sponge	Item code AHM0216 only: PVA (Polyvinyl Alcohol)			
	String	Item code AHM0216 only: 100% nylon silk thread			
Color	Gauze	White to light yellow			
	Sponge	Item code AHM0216 only: White			
	String	Item code AHM0216 only: Black			
Dimensions	Item Code		Gauze Size	Width (cm)	Length (cm)
	AHMDG20		0.5" x 0.5"	1,3	1,3
	AHM0022		2" x 2"	5,1	5,1
	AHM0024		2" x 4"	5,1	10,2
	AHM0044, AHM0444		4" x 4"	10,2	10,2
	Item Code	Sponge Size	Width (cm)	Length (cm)	Height (cm)
	AHM0216	2.25" x 1" x 0.6"	5,7	2,5	1,5
	Item Code	String Length (cm)			
	AHM0216	12,5			
Packaging	Item	Dressing or Gauze/Pouch	Pouches/ Box	Boxes/ Case	Dressings or Gauzes/Case
	AHM0022, AHM0024	1	12	80	960
	AHM0044		10		800
	AHM0216			10	100
	AHM0444	4	30	24	2280
	Item	Gauze/ Blister Pack	Blister Packs/ Box	Boxes/ Case	Gauzes/ Case
	AHMDG20	10	2	200	4000



TECHNICAL DATA SHEET
Trauma Gauze & Nasal Dressing
Items # AHM0022 to AHMDG20

Item Code	Description	Color	Specific Attributes
AHM0022	AMD-Ritmed Trauma Gauze, 5 x 5 cm (2" x 2"), 1 Gauze/Pouch, 12 Pouches/Box, 80 Boxes/Case, 960 Gauzes/Case	White	Sterile
AHM0024	AMD-Ritmed Trauma Gauze, 5 x 10 cm (2" x 4"), 1 Gauze/Pouch, 12 Pouches/Box, 80 Boxes/Case, 960 Gauzes/Case	White	Sterile
AHM0044	AMD-Ritmed Trauma Gauze, 10 x 10 cm (4" x 4"), 1 Gauze/Pouch, 10 Pouches/Box, 80 Boxes/Case, 800 Gauzes/Case	White	Sterile
AHM0216	AMD-Ritmed Nasal Dressing, 5,7 x 2,5 x 1,5 cm (2.25" x 1" x 0.6"), 1 Dressing/Pouch, 10 Pouches/Box, 10 Boxes/Case, 100 Dressings/Case	White	Sterile
AHM0444	AMD-Ritmed Trauma Gauze, 10 x 10 cm (4" x 4"), 4 Gauze/Pouch, 30 Pouches/Box, 24 Boxes/Case, 2880 Gauzes/Case	White	Sterile
AHMDG20	AMD-Ritmed Trauma Gauze, 1,27 x 1,27 cm (0.5" x 0.5"), 10 Gauzes/Blister Pack, 2 Blister Packs/Box, 200 Boxes/Case, 4000 Gauzes/Case	White	Sterile

Regulatory:

Not for sale in Canada (Class 2)
FDA Unclassified, Product Code FRO

Storage Conditions:

Store in a cool dry place

Country of Origin:

United States except item code AHM0216, China