

## Hemostyp Trauma Gauze

**SAFE, EFFECTIVE, AND EASY TO USE**

Hemostyp Trauma Gauze is a safe and cost effective tool used for the control of bleeding from the skin and other surface wounds where temporary control of bleeding is required.



### SAFE

Made from regenerated cotton cellulose  
Contains no animal proteins

### EFFECTIVE

When contacting blood and exudates, the product expands and forms into an adhesive gel  
Unique adherence and mimicking of the conformation of the wound site

### EASY TO USE

Easy application. Open the sterile pouch, fold or cut to best cover the bleeding surface, and apply with slight pressure  
Available in many convenient forms and sizes  
No prescription is required



<b>AHM0022</b>	2"x2", 1ply - trauma gauze	12 units per box
<b>AHM0024</b>	2"x4", 1ply - trauma gauze	12 units per box
<b>AHM0044</b>	4"x4", 1ply - trauma gauze	10 units per box
<b>AHM0216</b>	2.25"x1"x.6", 1ply - nasal dressing	10 units per box
<b>AHM0444</b>	4"x4", 1ply - trauma gauze	30 pouches of 4 units per box
<b>AHMDG20</b>	0.5"x0.5", 4ply - trauma gauze	20 units per box
<b>4900</b>	0.75"x0.75", 4ply - dental gauze	20 units per box