

MEDICOM® Technical Data Sheet

Fluid Impervious Gowns



Item Codes

Product Code	Description
8022	Regular - White
8024	X-Large - White

Intended Use


Non-surgical gown with waist ties, neck ties and knitted cuffs intended to protect health care personnel and patients from the transfer of microorganisms, body fluids (fluid repellent) and particulate material.

Packaging

10 gowns per bag
5 bags per case, 50 gowns per case

Size:

Area	Regular (cm)	Extra-Large (cm)
Gown Width (A)	71	78
Gown Length (B)	110	115
Sleeve Length (C)	58	61
Arm Width (D)	23	26
Arm Height (E)	28	32
Belt Length (H)	170	170



Color

White

Storage Conditions

Keep in a cool dry place

Material

Made of polypropylene with polyethylene coating.
Latex free.

Usage

Disposable, Single-Use only

Sterility

Non-Sterile

Country of Origin

China