

TECHNICAL DATA SHEET

Dissector Sponge Items # 53511 to 83560

Intended Use

Cotton dissector sponges which are used to separate or open vessels.

Characteristic		Specification					
Description		Gauze assembled in a spherical shape or woven tape rolled in a cylindrical shape with an X-Ray element. Held closed by stitches.					
Material		Gauze or Woven Tape	100% Cotton. Latex free.				
		X-Ray Element	Barium Sulfate (BaSO ₄) & PVC (Base material)				
Colour		Gauze or Woven Tape	White				
		X-Ray Element	Blue				
X-Ray Element		Position of X-Ray Element	Cylindrical	53511, 53513, 53531, 53535, 73512, 73514, 73518, 83516, 83522, 83532		Apparent, as a ring in the middle of the length	
			Spherical	53541, 53543, 85560		Apparent, circling the ball	
String		Item Code				String	
		53541, 53543, 83522, 83560				Single	
		53511, 53513, 53531, 53535, 73512, 73514, 73518, 83516, 83532				Non Strung	
Dimension		Size	Item Code	Shape	Diameter (mm)	Length (mm)	
		1/8" x 9/16"	53513, 73514	Cylindrical	5	15	
		1/4" x 9/16"	53511, 73512		8		
		1/4" x 3/4"	53535			25	
		1/4" x 1"	83516, 83522			19	
		1/2" x 1/2"	53531, 73518, 83532		13	13	
		1/4"	53541, 83560	Spherical	6	Not Applicable	
		1/2"	53543		13	Not Applicable	
Packaging	Non-Sterile Products	Item Code		Sponges per Glassine		Glassines per Case	Sponges per Case
		53531, 53535		5		100	500
		Item Code		Sponges per Count Card	Count Card per Glassine	Glassines per Case	Sponges per Case
		53541, 53543		5	1	100	500
		Item Code		Sponges per Holder	Holder per Glassine	Glassines per Case	Sponges per Case
		53511, 53513		5	1	100	500


TECHNICAL DATA SHEET

Dissector Sponge Items # 53511 to 83560

Characteristic		Specification					
Packaging	Sterile Products	Item Code	Sponges per Glassine	Glassine per Pouch	Pouches per Case	Sponges per Case	
		73518, 83516, 83532	5	1	20	100	
		Item Code	Sponges per Holder	Holder per Glassine	Glassine per Pouch	Pouches per Case	Sponges per Case
		73512, 73514	5	1	1	20	100
		Item Code	Sponges per Count Card	Count Card per Glassine	Glassine per Pouch	Pouches per Case	Sponges per Case
		83522	5	1	1	20	100
		83560				40	200

Visual

Cylindrical (53513)



Spherical (53543)



Item code: 83560 Spherical 'Rosebud'



Product Code	Description	Color	Specific Attributes
53511	AMD-Ritmed Dissector Sponge, $\frac{9}{16}$ " x $\frac{1}{4}$ ", 5 Sponges/Holder, 1 Holder/Glassine, 100 Glassines/Case, 500 Sponges/Case	White	Non-Sterile, Cylindrical, Non-Strung, X-ray, In Holder
53513	AMD-Ritmed Dissector Sponge, $\frac{9}{16}$ " x $\frac{1}{8}$ ", 5 Sponges/Holder, 1 Holder/Glassine, 100 Glassines/Case, 500 Sponges/Case	White	Non-Sterile, Cylindrical, Non-Strung, X-Ray, In Holder
53531	AMD-Ritmed Dissector Sponge, $\frac{1}{2}$ " x $\frac{1}{2}$ ", 5 Sponges/Glassine, 100 Glassines/Case, 500 Sponges/Case	White	Non-Sterile, Cylindrical, Non-Strung, X-Ray
53535	AMD-Ritmed Dissector Sponge, $\frac{1}{4}$ " x $\frac{3}{4}$ ", 5 Sponges/Glassine, 100 Glassines/Case, 500 Sponges/Case	White	Non-Sterile, Cylindrical, Non-Strung, X-Ray
53541	AMD-Ritmed Dissector Sponge, $\frac{1}{4}$ ", 5 Sponges/Count Card, 1 Count Card/Glassine, 100 Glassines/Case, 500 Sponges/Case	White	Non-Sterile, Spherical, Strung, X-Ray



TECHNICAL DATA SHEET

Dissector Sponge Items # 53511 to 83560

Product Code	Description	Color	Specific Attributes
53543	AMD-Ritmed Dissector Sponge, $\frac{1}{2}$ ", 5 Sponges/Count Card, 1 Count Card/Glassine, 100 Glassines/Case, 500 Sponges/Case	White	Non-Sterile, Spherical, Strung, X-Ray
73512	AMD-Ritmed Dissector Sponge, $\frac{9}{16}$ " x $\frac{1}{4}$ ", 5 Sponges/Holder, 1 Holder/Glassine, 1 Glassine/ Pouch, 20 Pouches/Case, 100 Sponges/Case	White	Sterile, Cylindrical, Non-Strung, X-Ray, In Holder
73514	AMD-Ritmed Dissector Sponge, $\frac{9}{16}$ " x $\frac{1}{8}$ ", 5 Sponges/Holder, 1 Holder/Glassine, 1 Glassine/ Pouch, 20 Pouches/Case, 100 Sponges/Case	White	Sterile, Cylindrical, Non-Strung, X-Ray, In Holder
73518	AMD-Ritmed Dissector Sponge, $\frac{1}{2}$ " x $\frac{1}{2}$ ", 5 Sponges/Glassine, 1 Glassine/Pouch, 20 Pouches/Case, 100 Sponges/Case	White	Sterile, Cylindrical, Non-Strung, X-Ray
83516	AMD-Ritmed Dissector Sponge, 1" x $\frac{1}{4}$ ", 5 Sponges/Glassine, 1 Glassine/Pouch, 20 Pouches/Case, 100 Sponges/Case	White	Sterile, Cylindrical, Non-Strung, X-Ray
83522	AMD-Ritmed Dissector Sponge, 1" x $\frac{1}{4}$ ", 5 Sponges/Glassine, 1 Glassine/Pouch, 20 Pouches/Case, 100 Sponges/Case	White	Sterile, Cylindrical, Strung, X-Ray
83532	AMD-Ritmed Dissector Sponge, $\frac{1}{2}$ " x $\frac{1}{2}$ ", 5 Sponges/Count Card, 1 Count Card/Glassine, 1 Glassine/Pouch, 20 Pouches/Case, 100 Sponges/Case	White	Sterile, Cylindrical, Non-Strung, X-Ray
83560	AMD-Ritmed Dissector Sponge Rosebud, $\frac{1}{4}$ ", 5 Sponges/Count Card, 1 Count Card/Glassine, 1 Glassine/Pouch, 40 Pouches/Case, 200 Sponges/Case	White	Sterile, Spherical, Strung, X-Ray

Regulatory

FDA Class 1, FDA Product Code GDY

Health Canada Class 2, Canadian Medical Device Number 5795

Storage Conditions

Store under ambient temperature and dry conditions

