

TECHNICAL DATA SHEET

Laboratory Coat Items # A8031 to A8051

Intended Use

Disposable SMS Laboratory Coat designed to protect health professionals, lab technicians and other users from fluids and cross contamination.

Characteristic	Specification		
Description	Non-woven gown with snap button front closure, knitted collar and cuffs. No pockets.		
Material	Spunbond Meltblown Spunbond (SMS) Polypropylene. Latex-Free		
Color	Item Code		Color
	A8031, A8032, A8033, A8034, A8035, A8036, A8038, A8039, A8040		Blue
	A8041, A8042, A8043, A8044, A8045, A8046, A8048, A8049, A8050, A8051		White
Dimension (Refer to Table 1)	Item Code		Size
	A8031, A8041		Small
	A8032, A8042		Medium
	A8033, A8043		Large
	A8034, A8044		XL
	A8035, A8045		XXL
	A8036, A8046		3XL
	A8038, A8048		4XL
	A8039, A8049		7XL
	A8040, A8050		9XL
	A8051		12XL
Packaging	Gowns / Bag	Bags / Case	Gowns / Case
	10	5	50

Table 1 – Gown Measurement

Dimension										
Measurements	Raglan Sleeve Style						Straight Sleeve Style			
	Small Size (cm)	Medium Size (cm)	Large Size (cm)	XL Size (cm)	XXL Size (cm)	3XL Size (cm)	4XL Size (cm)	7XL Size (cm)	9XL Size (cm)	12XL Size (cm)
(A) Underarm Width	62,5	64,5	66,5	69,0	72,0	73,5	77,0	86,0	95,0	101,0
(B) Lower Sleeve Length	48,0	49,0	51,0	52,0	53,0	55,0	57,5	65,0	70,0	68,0
(C) Cuff Length	5,0	5,0	5,0	5,0	5,0	5,0	5,0	5,0	5,0	5,0
(D) Shoulder Width	7,0	7,0	7,0	7,0	7,0	7,0	7,0	7,0	7,0	7,0
(E) Neck Width (stretched)	19,0	19,0	21,0	21,0	23,0	23,0	24,0	25,0	26,0	27,0
(F) Open Gown Width	134	139	144	149	154	158	164	182	198	212
(G) Gown Length	93,5	95,0	97,5	103,0	107,5	112,0	114,5	121,5	146,0	146,0
(H) Shoulder Height*	31,0	32,0	32,0	33,0	35,0	35,0	31,0	34,0	34,0	34,0
(I) Back Slit	30,5	30,5	30,5	30,5	30,5	30,5	30,5	30,5	30,5	30,5
(J) Sleeve Width at Wrist (stretched)	15,0	15,0	15,0	15,0	15,0	15,0	15,0	15,0	15,0	15,0
(K) Cuff Width	5,0	5,0	5,0	5,0	5,0	5,0	5,0	5,0	5,0	5,0
- Number of Buttons	5	5	5	5	5	6	6	6	6	6

TECHNICAL DATA SHEET

Laboratory Coat Items # A8031 to A8051

Figure 1 – Schematic representation of “raglan sleeve” gown (back view)





Figure 2 – Schematic representation of “straight sleeve” gown (back view)



TECHNICAL DATA SHEET

Laboratory Coat Items # A8031 to A8051

Item Code	Description	Color	Specific Attributes
A8031	AMD-Ritmed Laboratory Coat Small Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	Blue	Non-Sterile, Back Slit, Raglan Sleeve
A8032	AMD-Ritmed Laboratory Coat Medium Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	Blue	Non-Sterile, Back Slit, Raglan Sleeve
A8033	AMD-Ritmed Laboratory Coat Large Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	Blue	Non-Sterile, Back Slit, Raglan Sleeve
A8034	AMD-Ritmed Laboratory Coat XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	Blue	Non-Sterile, Back Slit, Raglan Sleeve
A8035	AMD-Ritmed Laboratory Coat XXL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	Bleu	Non-Sterile, Back Slit, Raglan Sleeve
A8036	AMD-Ritmed Laboratory Coat 3XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	Blue	Non-Sterile, Back Slit, Raglan Sleeve
A8038	AMD-Ritmed Laboratory Coat 4XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	Blue	Non-Sterile, Back Slit, Straight Sleeve
A8039	AMD-Ritmed Laboratory Coat 7XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	Blue	Non-Sterile, Back Slit, Straight Sleeve
A8040	AMD-Ritmed Laboratory Coat 9XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	Blue	Non-Sterile, Back Slit, Straight Sleeve
A8041	AMD-Ritmed Laboratory Coat Small Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	White	Non-Sterile, Back Slit, Raglan Sleeve
A8042	AMD-Ritmed Laboratory Coat Medium Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	White	Non-Sterile, Back Slit, Raglan Sleeve
A8043	AMD-Ritmed Laboratory Coat Large Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	White	Non-Sterile, Back Slit, Raglan Sleeve
A8044	AMD-Ritmed Laboratory Coat XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	White	Non-Sterile, Back Slit, Raglan Sleeve
A8045	AMD-Ritmed Laboratory Coat XXL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	White	Non-Sterile, Back Slit, Raglan Sleeve
A8046	AMD-Ritmed Laboratory Coat 3XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	White	Non-Sterile, Back Slit, Raglan Sleeve
A8048	AMD-Ritmed Laboratory Coat 4XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	White	Non-Sterile, Back Slit, Straight Sleeve
A8049	AMD-Ritmed Laboratory Coat 7XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	White	Non-Sterile, Back Slit, Straight Sleeve



TECHNICAL DATA SHEET

Laboratory Coat Items # A8031 to A8051

Item Code	Description	Color	Specific Attributes
A8050	AMD-Ritmed Laboratory Coat 9XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	White	Non-Sterile, Back Slit, Straight Sleeve
A8051	AMD-Ritmed Laboratory Coat 12XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	White	Non-Sterile, Back Slit, Straight Sleeve

Regulatory

FDA Class 1, FDA Product Code OEA
Health Canada Class 1

Storage Conditions

Store under ambient temperature and dry conditions

Country of origin

China