

TECHNICAL DATA SHEET

Laparotomy Sponges Items # 1501 to 1561K

Intended Use

Cotton gauze, with high absorbency and soft texture, used in operation room

Characteristic	Specification			
Description	Soft Cotton gauze with an X-ray and either no loop or one of the following attachments: loop, loop with plastic ring or loop with metal ring For Item Code 1501-ROLL: sponges are layered on top of one another and rolled. Apart from the regular X-Ray tape, an additional X-Ray element across the sponge is present			
Material	Fabric	100% Cotton (Latex Free)		
	X-Ray Element	Barium Sulfate (BaSO ₄) & PVC (Polyvinyl Chloride)		
	Tape Loop	100% Cotton (Latex Free)		
	Ring	Metal	Plastic	
		Spring Steel & Coated with Nickel on the surface	PE+ Barium sulfate	
Ply	4 ply			
Color	Fabric	White		
	X-Ray Element	Blue or Green		
	Tape Loop	Blue		
	Metal Ring	Silver		
	Plastic Ring	Grey		
Attachment	Type	Item Code		
	Tape Loop	1501, 1501-ROLL, 1502, 1503, 1503V, 1504, 1505, 1505-T, 1506, 1511, 1512, 1513, 1513-S, 1514, 1515, 1515-T, 1516, 1521, 1522, 1522-10, 1523, 1524, 1525, 1525-T, 1526, 1531, 1533, 1535, 1535-T, 1536, 1541-3, 1545, 1561, 1561K		
	Tape Loop with Metal Ring	1511-MR, 1521-MR		
	Tape Loop with Plastic Ring	1501-PR, 1511-PR, 1521-PR, 1531-PR, 1541-PR3, 1515R, 1515R-T		
	No Tape Loop	1501-NL, 1505-NL, 1511-NL, 1515NL, 1521-NL		
Dimension (Non-Washed Laparotomy Sponge)	Size	Item Code	Width (cm)	Length (cm)
	4" x 18"	1522, 1522-10, 1524, 1526	10.1	45.7
	8" x 36"	1536	20.3	91.4
	12" x 12"	1502, 1504, 1506,	30.4	30.4
	18" x 18"	1512, 1514, 1516	45.7	45.7

Laparotomy Sponges Items # 1501 to 1561K

Characteristic	Specification					
Dimension (Prewashed Laparotomy Sponge)	Size	Item Code	Non-Washed		Post Washed (stretched)	
			Width (cm)	Length (cm)	Width (cm)	Length (cm)
	4" x 18"	1521-MR, 1521-NL, 1521-PR, 1523, 1525, 1525-T	10.0	45.7	≥ 7.0	≥ 42.5
	8" x 36"	1531, 1531-PR, 1533, 1535, 1535-T	20.0	91.4	≥ 17.0	≥ 85.0
	12" x 12"	1501, 1501-NL, 1501-PR, 1501-ROLL, 1503, 1505, 1505-NL, 1505-T	30.4	30.4	≥ 26.0	≥ 27.0
	17" x 26"	1561, 1561K	43.1	66.0	≥ 42.5	≥ 61.5
	18" x 18"	1511, 1511-MR, 1511-NL, 1511-PR, 1513, 1513-S, 1515, 1515NL, 1515R, 1515R-T, 1515-T	45.7	45.7	≥ 40.0	≥ 42.5
	18" x 36"	1541-3, 1541-PR3, 1545		91.4	≥ 40.0	≥ 85.0
	1503V Diamond shaped with the tape loop located at the top of the design (left side)		Length (cm) (From tape loop corner to the opposite corner)			≥ 33.0
			Width (cm) (In the middle, perpendicular to length)			≥ 20.0
Ring	Dimension		Plastic		Metal	
	Outer Diameter (mm)		53		46	
	Inner Diameter (mm)		42		41	
Packaging	Item Code		Sponge per bundle	Bundle per bag	Bag per case	Sponge per case ⁽³⁾
	1541-PR3		3	36	2	216
	1541-3			40		240
	1503, 1504, 1513, 1513-S, 1514, 1523, 1524, 1533		5	20	2	200
	1511-MR, 1511-PR, 1531-PR, 1561, 1561K			40		400
	1511, 1511-NL, 1512, 1531			4	4	800
	1501-PR				2	700
	1501, 1501-NL, 1501-ROLL, 1502			4	4	1400
	1521-PR				2	1000
	1521-MR			1400		
	1521, 1521-NL, 1522			4	4	2800
	1522-10					

Laparotomy Sponges Items # 1501 to 1561K

Characteristic	Specification				
Packaging	Item Code	Sponge per bundle	Bundle per tray	Tray per case	Sponge per case ⁽³⁾
	1545	2	1	40	80
	1505, 1505-NL, 1505-T, 1506, 1515, 1515NL, 1515R, 1515R-T, 1515-T 1516, 1525, 1525-T, 1526, 1535, 1535-T, 1536	5			200
	Item Code	Sponge per pouch	Pouch per case		Sponge per case ⁽³⁾
	1503V	1	20		20

⁽¹⁾ N/A: Not Applicable

Item Code	Description	Color	Specific Attributes
1501	AMD-Ritmed Lap Sponges, 12"x12", 4 Ply, X-Ray, Tape Loop; Banded 5's, 70 pkg/paper bag, 4 bags/case, 1400 sponges/case	White	Prewashed, Steam Pretreatment, Non-Sterile
1501-NL	AMD-Ritmed Lap Sponges, 12"x12", 4 Ply, X-Ray, Banded 5's, No Loop, 70 pkg/paper bag, 4 bags/case, 1400 sponges/case	White	Prewashed, Steam Pretreatment, Non-Sterile
1501-PR	AMD-Ritmed Lap Sponges, 12"x12", 4 Ply, X-Ray, Tape Loop with Plastic Ring, Banded 5's, 70 pkg/paper bag, 2 bags/case, 700 sponges/case	White	Prewashed, Steam Pretreatment, Non-Sterile
1501-ROLL	AMD-Ritmed Lap Sponges, 12"x12", 4 Ply, X-Ray, Banded 5's, Rolled, Tape Loop, 70 pkg/paper bag, 4 bags/case, 1400 sponges/case	White	Prewashed, Steam Pretreatment, Non-Sterile
1502	AMD-Ritmed Lap Sponges, 12"x12", 4 Ply, X-Ray, Tape Loop, Banded 5's, 70 pkg/paper bag, 4 bags/case, 1400 sponges/case	White	Non-Washed, Steam Pretreatment, Non-Sterile
1503	AMD-Ritmed Lap Sponges, 12"x12", 4 Ply, X-Ray, Tape Loop; Banded 5's, 20 pkg/paper bag, 2 bags/case, 200 sponges/case	White	Prewashed, Steam Pretreatment, Non-Sterile
1503V	AMD-Ritmed Lap Sponges, V-Shape, X-Ray, Tape Loop; 1/glassine bag, 1 glassine bag/pouch, 20 sponges/case	White	Prewashed, Steam Pretreatment, Sterile
1504	AMD-Ritmed Lap Sponges, 12"x12", 4 Ply, X-Ray, Tape Loop; Banded 5's, 20 pkg/paper bag, 2 bags/case, 200 sponges/case	White	Non- Washed, Steam Pretreatment, Non-Sterile
1505	AMD-Ritmed Lap Sponges, 12"x12", 4 Ply, X-Ray, Tape Loop; Banded 5's, 40 trays/case, 200 sponges/case	White	Prewashed, Steam Pretreatment, Sterile
1505-NL	AMD-Ritmed Lap Sponges, 12"x12", 4 Ply, X-Ray, Banded 5's, No Loop, 40 trays/case, 200 sponges/case	White	Prewashed, Steam Pretreatment, Sterile
1505-T	AMD-Ritmed Lap Sponges, 12"x12", 4 Ply, X-Ray, Tape Loop, in Rigid Tray, Banded 5's, 40 trays/case, 200 sponges/case	White	Prewashed, Steam Pretreatment, Sterile
1506	AMD-Ritmed Lap Sponges, 12"x12", 4 Ply, X-Ray, Tape Loop, Banded 5's, 40 trays/case, 200 sponges/case	White	Non-Washed, Steam Pretreatment, Sterile
1511	AMD-Ritmed Lap Sponges, 18"x18", 4 Ply, X-Ray, Tape Loop; Banded 5's, 40 pkg/paper bag, 4 bags/case, 800 sponges/case	White	Prewashed, Steam Pretreatment, Non-Sterile

Laparotomy Sponges Items # 1501 to 1561K

Item Code	Description	Color	Specific Attributes
1511-MR	AMD-Ritmed Lap Sponges, 18"x18", 4 Ply, X-Ray, Tape Loop with Metal Ring, Banded 5's, 40 pkg/paper bag, 2 bags/case, 400 sponges/case	White	Prewashed, Steam Pretreatment, Non-Sterile
1511-NL	AMD-Ritmed Lap Sponges, 18"x18", 4 Ply, X-Ray, No Loop, Banded 5's, 40 pkg/paper bag, 4 bags/case, 800 sponges/case	White	Prewashed, Steam Pretreatment, Non-Sterile
1511-PR	AMD-Ritmed Lap Sponges, 18"x18", 4 Ply, X-Ray, Tape Loop with Plastic Ring, Banded 5's, 40 pkg/paper bag, 2 bags/case, 400 sponges/case	White	Prewashed, Steam Pretreatment, Non-Sterile
1512	AMD-Ritmed Lap Sponges, 18"x18", 4 Ply, X-Ray, Tape Loop, Banded 5's, 40 pkg/paper bag, 4 bags/case, 800 pcs/cs	White	Non-Washed, Steam Pretreatment, Non-Sterile
1513	AMD-Ritmed Lap Sponges, 18"x18", 4 Ply, X-Ray, Tape Loop, Banded 5's, 20 pkg/paper bag, 2 bags/case, 200 sponges/case	White	Prewashed, Steam Pretreatment, Non-Sterile
1513-S	AMD-Ritmed Lap Sponges, 18"x18", 4 Ply, X-Ray, Tape Loop, Banded 5's, 20 pkg/paper bag, 2 bags/case, 200 sponges/case	White	Prewashed, Steam Pretreatment, Gamma Pretreated
1514	AMD-Ritmed Lap Sponges, 18"x18", 4 Ply, X-Ray, Tape Loop, Banded 5's, 20 pkg/paper bag, 2 bags/case, 200 sponges/case	White	Non-Washed, Steam Pretreatment, Non-Sterile
1515	AMD-Ritmed Lap Sponges, 18"x18", 4 Ply, X-Ray, Tape Loop; Banded 5's, 40 trays/case, 200 sponges/case	White	Prewashed, Steam Pretreatment, Sterile
1515NL	AMD-Ritmed Lap Sponges, 18"x18", 4 Ply, X-Ray, No Loop, Banded 5's, 40 trays/case, 200 sponges/case	White	Prewashed, Steam Pretreatment, Sterile
1515R	AMD-Ritmed Lap Sponges, 18"x18", 4 Ply, X-Ray, Tape Loop with Plastic Ring, Banded 5's, 40 trays/case, 200 sponges/case	White	Prewashed, Steam Pretreatment, Sterile
1515R-T	AMD-Ritmed Lap Sponges, 18"x18", 4 Ply, X-Ray, Tape Loop with Plastic Ring, in Rigid Tray, Banded 5's, 40 trays/case, 200 sponges/case	White	Prewashed, Steam Pretreatment, Sterile
1515-T	AMD-Ritmed Lap Sponges, 18"x18", 4 Ply, X-Ray, Tape Loop, in Rigid Tray, Banded 5's, 40 trays/case, 200 sponges/case	White	Prewashed, Steam Pretreatment, Sterile
1516	AMD-Ritmed Lap Sponges, 18"x18", 4 Ply, X-Ray, Tape Loop; Banded 5's, 40 trays/case, 200 sponges/case	White	Not Washed, Steam Pretreatment, Sterile
1521	AMD-Ritmed Lap Sponges, 4"x18", 4 Ply, X-Ray, Tape Loop, Banded 5's, 140 pkg/paper bag, 4 bags/case, 2800 sponges/case	White	Prewashed, Steam Pretreatment, Non-Sterile
1521-MR	AMD-Ritmed Lap Sponges, 4"x18", 4 Ply, X-Ray, Tape Loop with Metal Ring, Banded 5's, 140 pkg/paper bag, 2 bags/case, 1400 sponges/case	White	Prewashed, Steam Pretreatment, Non Sterile
1521-NL	AMD-Ritmed Lap Sponges, 4"x18", 4 Ply, X-Ray, No Loop, Banded 5's, 140 pkg/paper bag, 4 bags/case, 2800 sponges/case	White	Prewashed, Steam Pretreatment, Non-Sterile
1521-PR	AMD-Ritmed Lap Sponges, 4"x18", 4 Ply, X-Ray, Tape Loop with Plastic Ring, Banded 5's, 100 pkg/paper bag, 2 bags/case, 1000 sponges/case	White	Prewashed, Steam Pretreatment, Non-Sterile
1522	AMD-Ritmed Lap Sponges, 4"x18", 4 Ply, X-Ray, Tape Loop, Banded 5's, 140 pkg/paper bag, 4 bags/case, 2800 sponges/case	White	Non-Washed, Steam Pretreatment, Non-Sterile
1522-10	AMD-Ritmed Lap Sponges, 4"x18", 4 Ply, X-Ray, Tape Loop, Banded 10's, 70 pkg/paper bag, 4 bags/case, 2800 sponges/case	White	Non-Washed, Steam Pretreatment, Non Sterile
1523	AMD-Ritmed Lap Sponges, 4"x18", 4 Ply, X-Ray, Tape Loop; Banded 5's, 20 pkg/paper bag, 2 bags/case, 200 sponges/case	White	Prewashed, Steam Pretreatment, Non-Sterile

Laparotomy Sponges Items # 1501 to 1561K

Item Code	Description	Color	Specific Attributes
1524	AMD-Ritmed Lap Sponges, 4"x18", 4 Ply, X-Ray, Tape Loop; Banded 5's, 20 pkg/paper bag, 2 bags/case, 200 sponges/case	White	Non-Washed, Steam Pretreatment, Non-Sterile
1525	AMD-Ritmed Lap Sponges, 4"x18", 4 Ply, X-Ray, Tape Loop, Banded 5's, 40 trays/case, 200 sponges/case	White	Prewashed, Steam Pretreatment, Sterile
1525-T	AMD-Ritmed Lap Sponges, 4"x18", 4 Ply, X-Ray, Tape Loop, in Rigid Tray, Banded 5's, 40 trays/case, 200 sponges/case	White	Prewashed, Steam Pretreatment, Sterile
1526	AMD-Ritmed Lap Sponges, 4"x18", 4 Ply, X-Ray, Tape Loop; Banded 5's, 40 trays/case, 200 sponges/case	White	Not Washed, Steam Pretreatment, Sterile
1531	AMD-Ritmed Lap Sponges, 8"x36", 4 Ply, X-Ray, Tape Loop, Banded 5's, 40 pkg/paper bag, 4 bags/case, 800 sponges/case	White	Prewashed, Steam Pretreatment, Non-Sterile
1531-PR	AMD-Ritmed Lap Sponges, 8"x36", 4 Ply, X-Ray, Tape Loop with Plastic Ring, Banded 5's, 40 pkg/paper bag, 2 bags/case, 400 sponges/case	White	Prewashed, Steam Pretreatment, Non-Sterile
1533	AMD-Ritmed Lap Sponges, 8"x36", 4 Ply, X-Ray, Tape Loop; Banded 5's, 20 pkg/paper bag, 2 bags/case, 200 sponges/case	White	Prewashed, Steam Pretreatment, Non-Sterile
1535	AMD-Ritmed Lap Sponges, 8"x36", 4 Ply, X-Ray, Tape Loop; Banded 5's, 40 trays/case, 200 sponges/case	White	Prewashed, Steam Pretreatment, Sterile
1535-T	AMD-Ritmed Lap Sponges, 8"x36", 4 Ply, X-Ray, Tape Loop; in Rigid Tray, Banded 5's, 40 trays/case, 200 sponges/case	White	Prewashed, Steam Pretreatment, Sterile
1536	AMD-Ritmed Lap Sponges, 8"x36", 4 Ply, X-Ray, Tape Loop; Banded 5's, 40 trays/case, 200 sponges/case	White	Non-Washed, Steam Pretreatment, Sterile
1541-3	AMD-Ritmed Lap Sponges, 18"x36", 4 Ply, X-Ray, Tape Loop, Banded 3's, 40 pkg/paper bag, 2 bags/case, 240 sponges/case	White	Prewashed, Steam Pretreatment, Non-Sterile
1541-PR3	AMD-Ritmed Lap Sponges, 18"x36", 4 Ply, X-Ray, Tape Loop with Plastic Ring, Banded 3's, 36 pkg/paper bag, 2 bags/case, 216 sponges/case	White	Prewashed, Steam Pretreatment, Non-Sterile
1545	AMD-Ritmed Lap Sponges, 18"x36", 4 Ply, X-Ray, Tape Loop; Banded 2's, 40 trays/case, 80 sponges/case	White	Prewashed, Steam Pretreatment, Sterile
1561K	AMD-Ritmed Lap Sponges, 17"x26", Fluff Gauze (Kerlix), X-Ray, Tape Loop, Banded 5's, 40 pkg/paper bag, 2 bags/case, 400 sponges/case	White	Non-Washed, Steam Pretreatment, Non-Sterile
1561	AMD-Ritmed Lap Sponges, 17"x26", 4 Ply, X-Ray, Tape Loop, Banded 5's, 40 pkg/paper bag, 2 bags/case, 400 sponges/case	White	Prewashed, Steam Pretreatment, Non-Sterile

Regulatory

FDA class 1 Product code GDY

Health Canada Class 2

Storage Conditions

Store under ambient temperature and dry conditions.

