


AMD-RITMED INC.			<div>Company Confidential</div> <div>Information contained in this document is the sole property of AMD-Ritmed Inc. Any unauthorized distribution or misuse is unlawful and may result in legal action(s).</div>
Drawing Title	3 ⁵ / ₈ ” Flat Wooden Spoon		
Technical Data Sheet			
Marketing: TDS 7.16-2 (PS26) Date: 29-May-2015	Unit of Measure: metric	Page 1 of 1	
Intended Use: Medical			
Regulatory: FDA Class 1, FDA Product Code ILC			



Packaging and Dimension								
Product Code	Pack (outer to inner most container)	Labeling		Length (L)	Width Bowl (W1)	Width Waist (W2)	Handle (H)	Thickness (T)
		Case	Inner Box					
51510	10/1000	Label	Printed	92.075 mm	23.00 mm	15.00 mm	16.891 mm	1.473 mm

Material	
Wooden Spoon	Paper (White) Birch, typically tan to light brown in color.
Inner Box	Non-Latex, Cardboard
Cases	Non-Latex, Corrugated

Specification	
Description	Flat piece of wood with rounded ends. One end larger than the other one. Unwrapped.
Color	Typically tan to light brown in color.