



TECHNICAL DATA SHEET

Non Adherent Pad Items # A1923 to A1938

Intended Use

Used for effective protection of cuts, scrapes and burns. Pads are easy to remove and don't stick to wound. Provide cushion protection.

Characteristic	Specification				
Description	Piece of fabric covered by a perforated film. Soft with smooth surface on both sides.				
Material Composition	Inside web: 80% cotton, 20% polyester Outer web: PET perforated film				
Color	White				
Dimension	Size	Item Code	Width (cm)	Length (cm)	
	2" x 3"	A1923, A1923-NS	5,1	7,6	
	3.5" x 8"	A1934	7,6	10,2	
	3" x 4"	A1938	8,9	20,3	
Packaging	Item Code	Pad/ Pouch	Pouches/ Box	Boxes/ Case	Pads/ Case
	A1923, A1934	1	100	12	1200
	A1938		50		600
	Item Code	Pads/Bag		Bags/Case	Pads/Case
	A1923-NS	100		72	7 200

Item Code	Description	Color	Specific Attributes
A1923	AMD-Ritmed Non Adherent Pads, 2"x 3", 1 Pad/Pouch, 100 Pouches/Box, 12 Boxes/Case, 1200 Pads/Case	White	Sterile
A1923-NS	AMD-Ritmed Non Adherent Pads, 2"x 3", Bulk, 100 Pads/ Bag, 72 Bags/Case, 7200 Pads/Case	White	Bulk, Non-Sterile
A1934	AMD-Ritmed Non Adherent Pads, 3"x 4", 1 Pad/Pouch, 100 Pouches/Box, 12 Boxes/Case, 1200 Pads/Case	White	Sterile
A1938	AMD-Ritmed Non Adherent Pads, 3 ½" x 8", 1 Pad/Pouch, 50 Pouches/Box, 12 Boxes/Case, 600 Pads/Case	White	Sterile

Regulatory

FDA Class 1, FDA Product Code NAC
Health Canada Class 1

Storage Conditions

Store under ambient temperature and dry conditions

Country of Origin

China