

TECHNICAL DATA SHEET

O.R. Towel Items # 5001-B to 5400-W

Intended Use

Cotton muslin Operation Room towels used to absorb biological fluids.

Characteristic	Specification			
Description	Flat muslin in the shape of a rectangle. Pre-washed items have a wrinkled appearance.			
Material Composition	Fabric	100% cotton, latex free		
	X-Ray Element	Barium Sulfate (BaSO ₄) & PVC (Base Material)		
X-Ray Element	For items which code ends with 'X': Located in the fold of the towel			
Dimension (pre-washed)	Size	Width (cm)	Length (cm)	
	17" x 25"	43,2	3,5	
Color	Item Code			Color
	5004-ND, 5400-ND			Non dyed
	5002-W, 5002-WX, 5004-W, 5004-WX, 5005-W, 5006-W, 5006-WX, 5400-W			White
	5001-B, 5002-B, 5002-BX, 5004-B, 5004-BX, 5006-B, 5008-B, 5010-B, 5100-B, 5100-BX, 5200-B, 5400-B			Blue
	5002-GX, 5004-G, 5006-G, 5006-GX, 5100-G			Green
Packaging	For sterile product only: CSR blue paper			
	Item Code	Towels/Pouch	Pouches/Case	Towel/Case
	5001-B	1	80	80
	5002-B, 5002-BX, 5002-GX, 5002-W, 5002-WX	2	40	
	5004-B, 5004-BX, 5004-G, 5004-ND, 5004-W, 5004-WX	4	20	
	5005-W	5	16	
	5006-B, 5006-G, 5006-GX, 5006-W, 5006-WX	6	12	72
	5008-B	8	10	80
	5010-B	10	8	
	Item Code	Towels/PE Bag	PE Bags/Case	Towels/Case
	5100-B, 5100-BX, 5100-G	100	1	100
	5200-B		2	200
	5400-B, 5400-ND, 5400-W		4	400



TECHNICAL DATA SHEET

O.R. Towel Items # 5001-B to 5400-W

Item Code	Description	Color	Specific Attributes
5001-B	AMD-Ritmed O.R. Towel, 1 Towel/Pouch, 80 Pouches/Case, 80 Towels/Case	Blue	Sterile
5002-B	AMD-Ritmed O.R. Towel, 2 Towels/Pouch, 40 Pouches/Case, 80 Towels/Case	Blue	Sterile
5002-BX	AMD-Ritmed O.R. Towel, 2 Towels/Pouch, 40 Pouches/Case, 80 Towels/Case	Blue	X-Ray Element, Sterile
5002-GX	AMD-Ritmed O.R. Towel, 2 Towels/Pouch, 40 Pouches/Case, 80 Towels/Case	Green	X-Ray Element, Sterile
5002-W	AMD-Ritmed O.R. Towel, 2 Towels/Pouch, 40 Pouches/Case, 80 Towels/Case	White	Sterile
5002-WX	AMD-Ritmed O.R. Towel, 2 Towels/Pouch, 40 Pouches/Case, 80 Towels/Case	White	X-Ray Element, Sterile
5004-B	AMD-Ritmed O.R. Towel, 4 Towels/Pouch, 20 Pouches/Case, 80 Towels/Case	Blue	Sterile
5004-BX	AMD-Ritmed O.R. Towel, 4 Towels/Pouch, 20 Pouches/Case, 80 Towels/Case	Blue	X-Ray Element, Sterile
5004-G	AMD-Ritmed O.R. Towel, 4 Towels/Pouch, 20 Pouches/Case, 80 Towels/Case	Green	Sterile
5004-ND	AMD-Ritmed O.R. Towel, 4 Towels/Pouch, 20 Pouches/Case, 80 Towels/Case	Non Dyed	Sterile
5004-W	AMD-Ritmed O.R. Towel, 4 Towels/Pouch, 20 Pouches/Case, 80 Towels/Case	White	Sterile
5004-WX	AMD-Ritmed O.R. Towel, 4 Towels/Pouch, 20 Pouches/Case, 80 Towels/Case	White	X-Ray Element, Sterile
5005-W	AMD-Ritmed O.R. Towel, 5 Towels/Pouch, 16 Pouches/Case, 80 Towels/Case	White	Sterile
5006-B	AMD-Ritmed O.R. Towel, 6 Towels/Pouch, 12 Pouches/Case, 72 Towels/Case	Blue	Sterile
5006-G	AMD-Ritmed O.R. Towel, 6 Towels/Pouch, 12 Pouches/Case, 72 Towels/Case	Green	Sterile
5006-GX	AMD-Ritmed O.R. Towel, 6 Towels/Pouch, 12 Pouches/Case, 72 Towels/Case	Green	X-Ray Element, Sterile
5006-W	AMD-Ritmed O.R. Towel, 6 Towels/Pouch, 12 Pouches/Case, 72 Towels/Case	White	Sterile
5006-WX	AMD-Ritmed O.R. Towel, 6 Towels/Pouch, 12 Pouches/Case, 72 Towels/Case	White	X-Ray Element, Sterile
5008-B	AMD-Ritmed O.R. Towel, 8 Towels/Pouch, 10 Pouches/Case, 80 Towels/Case	Blue	Sterile
5010-B	AMD-Ritmed O.R. Towel, 10 Towels/Pouch, 8 Pouches/Case, 80 Towels/Case	Blue	Sterile
5100-B	AMD-Ritmed O.R. Towel, 100 Towels/PE Bag, 1 PE Bag/Case, 100 Towels/Case	Blue	Non-Sterile
5100-BX	AMD-Ritmed O.R. Towel, 100 Towels/PE Bag, 1 PE Bag/Case, 100 Towels/Case	Blue	X-Ray Element, Non-Sterile
5100-G	AMD-Ritmed O.R. Towel, 100 Towels/PE Bag, 1 PE Bag/Case, 100 Towels/Case	Green	Non-Sterile
5200-B	AMD-Ritmed O.R. Towel, 100 Towels/PE Bag, 2 PE Bags/Case, 200 Towels/Case	Blue	Non-Sterile
5400-B	AMD-Ritmed O.R. Towel, 100 Towels/PE Bag, 4 PE Bags/Case, 400 Towels/Case	Blue	Non-Sterile
5400-ND	AMD-Ritmed O.R. Towel, 100 Towels/PE Bag, 4 PE Bags/Case, 400 Towels/Case	Non Dyed	Non-Sterile
5400-W	AMD-Ritmed O.R. Towel, 100 Towels/PE Bag, 4 PE Bags/Case, 400 Towels/Case	White	Non-Sterile



TECHNICAL DATA SHEET

O.R. Towel Items # 5001-B to 5400-W

Regulatory

FDA Class 2, Product Code KXX
Health Canada Class 1

Storage Condition

Store under ambient temperature and dry conditions

Country of Origin

China