

Safe+Mask® Cone Masks with Earloop, Non Sterile #2020 & #2021

Material characteristics	100% Polyester Latex Free
Weight	200 g/m ²
Elastic Band Composition	White Nylon Elastic (ultrasonic welding)
Nose-piece	Aluminum
Color	Available in blue, pink
Usage	Single Use Only
Appearance	Molded contour mask (“cone” shape) with nosepiece and stretchable headband.
Intended Use	The cone masks are intended to protect nose and mouth of the wearer during general medical procedures.



Dimensions			Packaging
Mask Body	Length	117 – 137 mm	50 masks / box 12 boxes / case 600 masks /case
	Width	104 – 117 mm	
	Height	40 – 60 mm	
Nose Piece	Length	77 – 81 mm	
Elastic Band	Length	285 – 350 mm	

Storage	Store under cool, clean, and dry conditions. Avoid excessive heat (over 40°C or 104°F)
----------------	--