



TECHNICAL DATA SHEET

Neurological Sponge, Single X-Ray, Single Strung Items # 55903C to 85999C

Intended Use

Sponges used to moisten or to absorb biological fluid during neurological surgery. May also be used to protect brain during incision.

Characteristic	Acceptable Limits									
Description	Soft non woven fabric									
Material Composition	Fabric	Polyester and viscose								
	String	Polyester								
	X-Ray Element	Barium Sulfate (BaSO ₄) & PVC (Base Material)								
String	Single strung									
X-Ray Element	Single X-Ray									
Color	Sponge		String		X-Ray Element					
	White		Green		Blue					
Dimension	Size	Item Code		Width (cm)	Length (cm)					
	1/4" x 1/4"	55991C, 85991C		0,7	0,7					
	1/4" x 1 1/4"	55903C			3,2					
	1/4" x 3"	55921C, 85936C			7,6					
	1/2" x 1/2"	55992C, 85992C		1,3	1,3					
	1/2" x 1"	55993C, 85993C			2,6					
	1/2" x 1 1/2"	55998C, 85998C			3,8					
	1/2" x 2"	55999C, 85999C			5,1					
	1/2" x 3"	55994C, 85994C			7,6					
	3/4" x 3"	55922C, 85938C		2,0	7,6					
Packaging	Item Code		Sponges/ Glassine		Glassines/ Case		Sponges/ Case ⁽⁴⁾			
	55903C, 55921C, 55922C, 55991C, 55992C, 55993C, 55994C, 55998C, 55999C		10		10		100			
	Item Code		Sponges/ Glassine		Glassine/ Pouch		Pouch/ Case		Sponges/ Case ⁽⁴⁾	
	85936C, 85938C, 85991C, 85992C, 85993C, 85994C, 85998C, 85999C		10		1		10		100	

Item Code	Description	Color	Specific Attributes
55903C	AMD-Ritmed Neurological Sponge, 1/4"x1 1/4", 10 Sponges/Glassine, 10 Glassines/Case, 100 Sponges/Case	White	Single X-Ray Element, Single Strung, Non-Sterile
55921C	AMD-Ritmed Neurological Sponge, 1/4"x3", 10 Sponges/Glassine, 10 Glassines/Case, 100 Sponges/Case	White	Single X-Ray Element, Single Strung, Count Card, Non-Sterile



TECHNICAL DATA SHEET

Neurological Sponge, Single X-Ray, Single Strung Items # 55903C to 85999C

Item Code	Description	Color	Specific Attributes
55922C	AMD-Ritmed Neurological Sponge, ¾"x3", 10 Sponges/Glassine, 10 Glassines/Case, 100 Sponges/Case	White	Single X-Ray Element, Single Strung, Count Card, Non-Sterile
55991C	AMD-Ritmed Neurological Sponge, ¼"x¼", 10 Sponges/Glassine, 10 Glassines/Case, 100 Sponges/Case	White	Single X-Ray Element, Single Strung, Count Card, Non-Sterile
55992C	AMD-Ritmed Neurological Sponge, ½"x½", 10 Sponges/Glassine, 10 Glassines/Case, 100 Sponges/Case	White	Single X-Ray Element, Single Strung, Count Card, Non-Sterile
55993C	AMD-Ritmed Neurological Sponge, ½"x1", 10 Sponges/Glassine, 10 Glassines/Case, 100 Sponges/Case	White	Single X-Ray Element, Single Strung, Count Card, Non-Sterile
55994C	AMD-Ritmed Neurological Sponge, ½"x3", 10 Sponges/Glassine, 10 Glassines/Case, 100 Sponges/Case	White	Single X-Ray Element, Single Strung, Count Card, Non-Sterile
55998C	AMD-Ritmed Neurological Sponge, ½"x1½", 10 Sponges/Glassine, 10 Glassines/Case, 100 Sponges/Case	White	Single X-Ray Element, Single Strung, Count Card, Non-Sterile
55999C	AMD-Ritmed Neurological Sponge, ½"x2", 10 Sponges/Glassine, 10 Glassines/Case, 100 Sponges/Case	White	Single X-Ray Element, Single Strung, Count Card, Non-Sterile
85936C	AMD-Ritmed Neurological Sponge, ¼"x3", 10 Sponges/Glassine, 1 Glassine/Pouch, 10 Pouches/Case, 100 Sponges/Case	White	Single X-Ray Element, Single Strung, Count Card, Sterile
85938C	AMD-Ritmed Neurological Sponge, ¾"x3", 10 Sponges/Glassine, 1 Glassine/Pouch, 10 Pouches/Case, 100 Sponges/Case	White	Single X-Ray Element, Single Strung, Count Card, Sterile
85991C	AMD-Ritmed Neurological Sponge, ¼"x¼", 10 Sponges/Glassine, 1 Glassine/Pouch, 10 Pouches/Case, 100 Sponges/Case	White	Single X-Ray Element, Single Strung, Count Card, Sterile
85992C	AMD-Ritmed Neurological Sponge, ½"x½", 10 Sponges/Glassine, 1 Glassine/Pouch, 10 Pouches/Case, 100 Sponges/Case	White	Single X-Ray Element, Single Strung, Count Card, Sterile
85993C	AMD-Ritmed Neurological Sponge, ½"x1", 10 Sponges/Glassine, 1 Glassine/Pouch, 10 Pouches/Case, 100 Sponges/Case	White	Single X-Ray Element, Single Strung, Count Card, Sterile



TECHNICAL DATA SHEET

Neurological Sponge, Single X-Ray, Single Strung Items # 55903C to 85999C

Item Code	Description	Color	Specific Attributes
85994C	AMD-Ritmed Neurological Sponge, ½"x3", 10 Sponges/Glassine, 1 Glassine/Pouch, 10 Pouches/Case, 100 Sponges/Case	White	Single X-Ray Element, Single Strung, Count Card, Sterile
85998C	AMD-Ritmed Neurological Sponge, ½"x1½", 10 Sponges/Glassine, 1 Glassine/Pouch, 10 Pouches/Case, 100 Sponges/Case	White	Single X-Ray Element, Single Strung, Count Card, Sterile
85999C	AMD-Ritmed Neurological Sponge, ½"x2", 10 Sponges/Glassine, 1 Glassine/Pouch, 10 Pouches/Case, 100 Sponges/Case	White	Single X-Ray Element, Single Strung, Count Card, Sterile

Regulatory

FDA Class 1, FDA Product Code GDY, Listing Number E637831

Health Canada Class 2, Canadian License Number 4348

Storage Conditions

Store under ambient temperature and dry conditions