

TECHNICAL DATA SHEET

Countway System Items # 2213 to 2214

Intended Use

Single use medical device which is used in operation rooms (O.R.) to handle soiled sponges and determine blood loss. This system has unique clear self-closing pouches that contains and confines soiled sponges and makes counting quick, simple and accurate. The covered pouches prevent odour and decreases the risk of airborne contamination protecting healthcare personnel and patient and maintaining better O.R. sanitation. Pockets are clear to enable visibility and facilitate counting. Single use


Characteristic	Acceptable Limits				
Description	2213, 2213-BL and 2213-NS: Countway system panel with 5 self-closing pockets (divided or not) on one side. Pockets are clear to enable visibility and facilitate counting 2214: Countway system inner bag which is made of clear plastic				
Material Composition	Polyethylene				
Color	Item Code		Color		
	2213, 2213-NS		Green		
	2213-BL		Blue		
	2214		Clear		
Dimension	Item Code	Part	Size	Width (cm)	Length (cm)
	2213, 2213-NS,	Panel	14" x 23 ¼"	35,5	59
		Panel Flap	14" x 3 ½"		9
		Pocket	6 ¾" x 4 ½"	17	11,5
		Pocket Flap	7" x 5"	18	12,5
	2213-BL	Panel	14" x 30 ¾"	35,5	78
		Panel Flap	14" x 3 ½"		9
		Pocket	6" x 6 ¾"	15	17
		Pocket Flap	6 ½" x 7"	16,5	18
	2214	Part	Open Bag Width (cm)		Height (cm)
		Bag	61		51
Packaging	Item Code		Bags or Panels / Box	Box / Case	Bags or Panels / Case
	2213, 2213-BL, 2213-NS, 2214		50	8	400



TECHNICAL DATA SHEET

Countway System Items # 2213 to 2214

Countway Diagram



Item Code	Description	Color	Specific Attributes
2213	AMD-Ritmed Countway System Panel, 50 Panels/Box, 8 Boxes/Case, 400 Panels/Case	Green	5 Pockets, Divider
2213-BL	AMD-Ritmed Countway System Panel, 50 Panels/Box, 8 Boxes/Case, 400 Panels/Case	Blue	5 pockets, No Divider
2213-NS	AMD-Ritmed Countway System Panel, 50 Panels/Box, 8 Boxes/Case, 400 Panels/Case	Green	5 pockets, No Divider
2214	AMD-Ritmed Countway System Inner Bag, 50 Panels/Box, 8 Boxes/Case, 400 Panels/Case	Clear	N/A

Regulatory

FDA Class 1, FDA Product Code LWH

Health Canada, Class 1

Storage Conditions

Store under ambient and dry conditions

Country of Origin

China except 2214, Canada