

TECHNICAL DATA SHEET

Gauze Sponge (All Gauze) Items # A2120 to A8800

Intended Use

May be used for one or more of the following purposes, cleansing wounds or skin, absorbing biological fluids or support organs and tissue during surgical procedures

Characteristic	Specification					
Description	Folded gauze					
Material Composition	100% cotton ⁽¹⁾ , Latex free					
Color	White					
Slit	For Tracheotomy Sponges only: A slit is present between the center of the Sponge and the folded edge. Slit edge is sewn inward					
Ply	Item Code			Folding ⁽²⁾	Ply	
	A2800, A3001, A3001-10, A3003, A3005, A3005-10, A3008, A3013, A3800, A4800, A5800, A8800			C	8	
	A2120, A2212, A3002, A3004, A3006, A3120, A3312, A4120, A4412				12	
	A3009, A3010, A4412T, A8120, A8121, A8125			E		
	A4121, A4121-T, A4125,			C	16	
	A3007, A3011, A3026SX5, A4161					
Dimension	Size	Item Code		Length (cm)	Width (cm)	
	2"x2"	A2120, A2212, A2800, A3001, A3001-10, A3002		4,5 – 5,5	4,5 – 5,5	
	3"x3"	A3003, A3004, A3120, A3312, A3800		7,0 – 8,0	7,0 – 8,0	
	3"x4"	A3013, A5800			9,5 – 10,5	
	4"x4"	A3005, A3005-10, A3006, A3007, A3010, A3011, A3026SX5, A4120, A4121, A4121-T, A4125, A4161, A4412, A4412T, A4800		9,5 – 10,5		
	8"x4"	A8120, A8121, A8125, A8800		19,2 – 21,2		
		A3008, A3009		19,0 – 20,0		
Packaging	Item Code		Sponges/Bag	Bags/Case	Sponges/Case	
	A3006, A3007, A3008, A3009, A3010, A3011, A3013		200	10	2000	
	A3003, A3004, A3005, A3005-10			20	4000	
	A3001, A3001-10			25	5000	
	A3002			40	8000	
	Item Code		Pads/Pouch	Pouches/Box	Boxes/Case	Pads/Case
	4412		1	100	12	2400
	A2212, A3312				24	2400

TECHNICAL DATA SHEET

Gauze Sponge (All Gauze) Items # A2120 to A8800


Characteristic	Specification				
Packaging	Item Code	Sponges/Pouch	Pouches/Tray	Trays/ Case	Sponges/Case
	A4412T	1	25	24	600
	A8120, A8800	2		20	1000
	A4120			24	1200
	A3120, A3800		40	30	2400
	A4800, A5800		50	12	1200
	A2120, A2800			30	3000
	A8125	5	10	20	1000
	A4125		20	12	1200
	A4121, A4161, A8121	10	10		1200
	Item Code	Sponges/Hard Plastic Tray	Hard Plastic Trays/Box	Boxes/Case	Sponges/Case
	A4121-T	10	64	2	1280
	Item Code	Sponges/Pouch	Pouches/Case		Sponges/Case
	A3026SX5	5	100		500

(1) The cotton yarn used for all woven gauze products are Kier bleached


(2) Refer to folding diagram

Folding Diagram


C Fold - 8 ply



C Fold - 12 ply




E Fold - 12 ply



TECHNICAL DATA SHEET

Gauze Sponge (All Gauze) Items # A2120 to A8800

C Fold -16 ply



Item Code	Description	Color	Specific Attribute
A2120	AMD-Ritmed Gauze Sponges, 2"x2", 12 Ply, 2 Sponges/Pouch, 50 Pouches/Tray, 30 Trays/Case, 3000 Sponges/Case	White	Sterile
A2212	AMD-Ritmed Gauze Pads, 2"x2", 12 Ply, 1 Sponge/Pouch, 100 Pouches/box, 24 Pouches/box, 2400 Sponges/Case	White	Sterile
A2800	AMD-Ritmed Gauze Sponges, 2"x2", 8 Ply, 2 Sponges/Pouch, 50 Pouches/box, 30 boxes/Case, 3000 Sponges/Case	White	Sterile
A3001	AMD-Ritmed Gauze Sponges, 2"x2", 8 Ply, 200 Sponges/Bag, 25 Bags/Case, 5000 Sponges/Case	White	Non-Sterile
A3001-10	AMD-Ritmed Gauze Sponges, 2"x2", 8 Ply, 200 Sponges/Bag, 25 Bags/Case, 5000 Sponges/Case	White	Non-Sterile, Indexed in 10's(*)
A3002	AMD-Ritmed Gauze Sponges, 2"x2", 12 Ply, 200 Sponges/Bag, 40 Bags/Case, 8000 Sponges/Case	White	Non-Sterile
A3003	AMD-Ritmed Gauze Sponges, 3"x3", 8 Ply, 200 Sponges/Bag, 20 Bags/Case, 4000 Sponges/Case	White	Non-Sterile
A3004	AMD-Ritmed Gauze Sponges, 3"x3", 12 Ply, 200 Sponges/Bag, 20 Bags/Case, 4000 Sponges/Case	White	Non-Sterile
A3005	AMD-Ritmed Gauze Sponges, 4"x4", 8 Ply, 200 Sponges/Bag, 20 Bags/Case, 4000 Sponges/Case	White	Non-Sterile
A3005-10	AMD-Ritmed Gauze Sponges, 4"x4", 8 Ply, 200 Sponges/Bag, 20 Bags/Case, 4000 Sponges/Case	White	Non-Sterile, Indexed in 10's(*)
A3006	AMD-Ritmed Gauze Sponges, 4"x4", 12 Ply, 200 Sponges/Bag, 10 Bags/Case, 2000 Sponges/Case	White	Non-Sterile
A3007	AMD-Ritmed Gauze Sponges, 4"x4", 16 Ply, 200 Sponges/Bag, 10 Bags/Case, 2000 Sponges/Case	White	Non-Sterile
A3008	AMD-Ritmed Gauze Sponges, 8"x4", 8 Ply, 200 Sponges/Bag, 10 Bags/Case, 2000 Sponges/Case	White	Non-Sterile
A3009	AMD-Ritmed Gauze Sponges, 8"x4", 12 Ply, 200 Sponges/Bag, 10 Bags/Case, 2000 Sponges/Case	White	Non-Sterile
A3010	AMD-Ritmed Gauze Sponges, 4"x4", 12 Ply, 200 Sponges/Bag, 10 Bags/Case, 2000 Sponges/Case	White	Tracheotomy Sponge, Non-Sterile
A3011	AMD-Ritmed Gauze Sponges, 4"x4", 16 Ply, 200 Sponges/Bag, 10 Bags/Case, 2000 Sponges/Case	White	Tracheotomy Sponge, Non-Sterile
A3013	AMD-Ritmed Gauze Sponges, 3"x4", 8 Ply, 200 Sponges/Bag, 10 Bags/Case, 2000 Sponges/Case	White	Non-Sterile
A3026SX5	AMD-Ritmed Gauze Sponges, 4"x4", 16 Ply, 5 Sponges/Pouch, 100 Pouches/Case, 500 Sponges/Case	White	Sterile

TECHNICAL DATA SHEET

Gauze Sponge (All Gauze) Items # A2120 to A8800

Item Code	Description	Color	Specific Attribute
A3120	AMD-Ritmed Gauze Sponges, 3"x3", 12 Ply, 2 Sponges/Pouch, 40 Pouches/Tray, 30 Trays/Case, 2400 Sponges/Case	White	Sterile
A3312	AMD-Ritmed Gauze Pads, 3"x3", 12 Ply, 1 Sponge/Pouch, 100 Pouches/Tray, 24 Trays/Case, 2400 Sponges/Case	White	Sterile
A3800	AMD-Ritmed Gauze Sponges, 3"x3", 8 Ply, 2 Sponges/Pouch, 40 Pouches/Tray, 30 Trays/Case, 2400 Sponges/Case	White	Sterile
A4120	AMD-Ritmed Gauze Sponges, 4"x4", 12 Ply, 2 Sponges/Pouch, 25 Pouches/Tray, 24 Trays/Case, 1200 Sponges/Case	White	Sterile
A4121	AMD-Ritmed Gauze Sponges, 4"x4", 12 Ply, 10 Sponges/Pouch, 10 Pouches/Tray, 12 Trays/Case, 1200 Sponges/Case	White	Sterile
A4121-T	AMD-Ritmed Gauze Sponges, 4"x4", 12 Ply, 10 Sponges/Tray, x 128 Trays/Case, 1280 Sponges/Case	White	Sterile, Hard Plastic Tray
A4125	AMD-Ritmed Gauze Sponges, 4"x4", 12 Ply, 5 Sponges/Pouch, 20 Pouches/Tray, 12 Trays/Case, 1200 Sponges/Case	White	Sterile
A4161	AMD-Ritmed Gauze Sponges, 4"x4", 16 Ply, 10 Sponges/Pouch, 10 Pouches/Tray, 12 Trays/Case, 1200 Sponges/Case	White	Sterile
A4412	AMD-Ritmed Gauze Pads, 4"x4", 12 Ply, 1 Sponge/Pouch, 100 Pouches/Tray, 12 Trays/Case, 1200 Sponges/Case	White	Sterile
A4412T	AMD-Ritmed Gauze Sponges, 4"x4", 12 Ply, 1 Sponge/Pouch, 25 Pouches/Tray, 24 Trays/Case, 600 Sponges/Case	White	Tracheotomy Sponge, Sterile
A4800	AMD-Ritmed Gauze Sponges, 4"x4", 8 Ply, 2 Sponges/Pouch, 50 Pouches/Tray, 12 Trays/Case, 1200 Sponges/Case	White	Sterile
A5800	AMD-Ritmed Gauze Sponges, 3"x4", 8 Ply, 2 Sponges/Pouch, 50 Pouches/Tray, 12 Trays/Case, 1200 Sponges/Case	White	Sterile
A8120	AMD-Ritmed Gauze Sponges, 8"x4", 12 Ply, 2 Sponges/Pouch, 25 Pouches/Tray, 20 Trays/Case, 1000 Sponges/Case	White	Sterile
A8121	AMD-Ritmed Gauze Sponges, 8"x4", 12 Ply, 10 Sponges/Pouch, 10 Pouches/Tray, 12 Trays/Case, 1200 Sponges/Case	White	Sterile
A8125	AMD-Ritmed Gauze Sponges, 8"x4", 12 Ply, 5 Sponges/Pouch, 10 Pouches/Tray, 20 Trays/Case, 1000 Sponges/Case	White	Sterile
A8800	AMD-Ritmed Gauze Sponges, 8"x4", 8 Ply, 2 Sponges/Pouch, 25 Pouches/Tray, 20 Trays/Case, 1000 Sponges/Case	White	Sterile

* Every 10th (-10) gauze is not aligned with the rest of the gauzes in the package in order to serve as an indicator

Regulatory

FDA Class 1, FDA Product Code NAC

Health Canada Class 1

Storage Conditions

Store under ambient temperature and dry conditions

Country of Origin

China except for item code A3026SX5, Canada