


AMD-RITMED INC.		
Title	Wood (Dowel) Applicator	
Technical Data Sheet		
Marketing: TDS 7.12-1 (PS4) Date: 29-May-2015	Unit of Measure: Inches	Page 1 of 1
Intended Use: Medical - General Hospital		
Regulatory: FDA Class 1, FDA Product Code JRQ Country of Origin: China		

Company Confidential

*Information contained in this document is the sole property of AMD-Ritmed Inc.
Any unauthorized distribution or misuse is unlawful and may result in legal action(s).*

Operation Management Approval:



Dry Dimension				
Item Code	Length (L)	Diameter (D)	Maximum Out of Round	Maximum Diameter Variation
55900	5.813" (147.6 mm)	0.085" (2.16 mm)	0.003" (0.08mm)	0.005" (0.13 mm)

Material	
Applicator	Paper (White) Birch
	Color: Tan with natural yellow or reddish, brown tints
Inner Boxes	Non Latex, Box
Cases	Non-Latex, Corrugated Box

Product Configuration Table					
Product	Sterility	** Packaging		* Labeling	
Item Code	Non Sterile	Inner Box (Clay Coat)	Pack (Outer to Inner most container)	Case (Outer Box)	Inner Box
55900	√	√	20/1000	L	P

* P = Printed Container	L = Labeled Container
**All wood applicators are packed in an outer corrugated case	