

Universal Drainage Bag

Versatility. Durability. Flexibility.



- Flexible silicone tubing facilitates trimming to desired length.
- Soft Velcro® strap cushions securement of bag system to patient.
- Twist valve facilitates quick and easy disposal of bag contents.
- Reflux valve prevents retrograde contamination.

Universal Drainage Bag

Contents:

1 ea. Universal Drainage bag
1 ea. Velcro® Strap (Individually Packaged)

- 600 ml PVC bag
(4 13/16" X 10 5/8")
with 100 ml graduations
- 20" long drainage tube
with male luer adapter at proximal
end



Reorder Number	Product Code	Units/case
622290	FCS-2000	20

Ordering Information:

For more information or to place an order,
contact your B. Braun Interventional Systems
representative or call **1-877-VENA CAV (836-2228)**

