



## Features and Benefits of the Heart Sync Product Line

### Cost Savings

On average, we can save an organization up to 15% on their annual defibrillation electrode spend.

### Standardization

We can standardize an organization to one electrode from one vendor no matter what make or model of defibrillator you have.

### Three Year Shelf Life

Extended shelf life reduces cost and waste. Our Adult/Child and Pediatric electrodes in low use areas will provide an additional 50% shelf life versus the rest of the market.

### Longest Lead Wire on the Market

60" lead wire. The longer lead wire gives the clinician more flexibility while positioning the defibrillator. This is also beneficial with larger patients.

### Color Keyed Pouches

Eliminates mistakes. Our labels are consistent with the different connector styles: Physio in red, Philips in green, Zoll in yellow, and the Anderson label is blue. Our Adult/Child pouch background is white and our Pediatric pouch background is light blue offering clinicians a visual cue between the two products.

### Peel Pouch Bags

Easier for clinicians to open. Pouch also has traditional tear notch at the bottom.

### User Instructions on Release Liner

User instructions and a check list is on the colored release liner. A colored liner does not become a slipping hazard for the staff while on the treatment floor.

### Smaller Pad Size

Increasing the conductor plate while reducing the overall size of the electrode. Our aggressive conductive hydrogel reduces the need for a large foam border to anchor the pad to the patient.

### Color Keyed Pads

Anterior electrode is blue foam and the posterior electrode is white foam. This helps train the staff to look for colors to aid in pad placement.

### Instructions on the Package

Pouch instructions designed only for the US market. Complete user instructions, cautions and warnings on the label, not an IFU in the box.