

Micro-Touch® E.P.® Gloves

Micro-Touch® E.P.® Gloves are made from natural rubber latex, feature extra long reinforced cuffs, and are thicker for extra protection, making them ideal for use in higher risk clinical applications. These powder-free gloves are hand specific, which allow for added comfort and fit.

Safety & Protection

Infection Control

Allergy Concerns. This examination glove contains natural rubber latex which may cause allergic reactions. Safe use of these gloves by latex-sensitized individuals has not been established.

Sterilization. This is a non-sterile examination glove.

Quality Benchmarks. Meets or exceeds ASTM examination glove standards. Manufactured within the quality guidelines of ISO 13485:2003, FDA-QSR, and CE Certification. FDA approved for handling chemotherapy drugs.¹

Cuff. Extended 12" straight, reinforced cuff.

Resistance. Latex remains the gold standard for hand barrier protection due to its strength, proven barrier protection, elasticity, fit, feel, and comfort. This glove is twice as thick as standard latex gloves for added protection in higher-risk clinical settings.

This glove is also tested for use with the following 18 chemotherapy drugs:

Bleomycin Sulfate, Carboplatin, Carmustine, Cisplatin, Cyclophosphamide, Cytarabine Hydrochloride, Daunorubicin Hydrochloride, Doxorubicin Hydrochloride, Etoposide, 5-Fluorouracil, Idarubicin, Ifosfamide, Mechlorethamine Hydrochloride, Melphalan, Mitomycin, Mitoxantrone, Paclitaxel, and ThioTEPA.

Skin Protection

Powder-Free. To help eliminate powder-related complications and sensitivity.

Comfort & Fit

Hand Fatigue

Stretchability (Modulus). Latex has excellent elasticity, allowing the glove to stretch easily and return to its original shape.

Shape of Glove/Former. Hand specific for excellent comfort and fit.

Tactile Sensitivity

Dependability. Because this glove is latex, wearers experience a more secure fit that provides greater comfort without bagginess.

Memory. Latex has excellent elasticity and memory.

Grip

Application Driven. The fully textured finish of the glove provides extra protection for wet or dry grip and against slippage.

Ordering Information

SIZE	Product Number
X-Small	9201
Small	9203
Medium	9205
Large	9207

1 Pair/Package

25 Pairs/Dispenser Box – 4 Boxes/Case



Caution: This product contains natural rubber latex, which may cause allergic reactions.

1. Specific chemical resistance data available upon request: 1-800-952-9916 (U.S.A. only). Gloves used for protection against chemotherapy drugs must be selected specifically for the type of chemicals used. Review material safety data sheets for the chemicals being used to determine the required level of protection.

MICRO-TOUCH®
E.P.® GLOVES



Ansell

Everything you touch...we touch.