



MICRO-TOUCH[®]



EXAMINATION AND GENERAL PURPOSE GLOVE CATALOG

Ansell

Everything you touch...we touch.



Table of Contents

Introduction	4-5
Latex-Free, Powder-Free Examination Gloves	6-8
Nitrile	6-7
Micro-Touch® NitraFree™	6
Micro-Touch® Nitrile	6
Micro-Touch® Smooth Nitrile	7
Micro-Touch® Nitrile E.P.®	7
Neoprene	8
Micro-Touch® Affinity®	8
Vinyl	8
Micro-Touch® Elite®	8
Sterile, Powder-Free Examination Gloves	9
Micro-Touch® Nitra-Tex® Sterile Pairs and Singles	9
Micro-Touch® Plus Sterile Singles	9
Latex, Powder-Free Examination Gloves	10-11
Micro-Touch® Plus	10
Micro-Touch® NextStep®	10
Micro-Touch® E.P.® Gloves	11
Latex, Powdered Examination Gloves	12
Micro-Touch® Powder	12
General Purpose Gloves	13-14
Housekeeping	13
Latex/Nitrile Utility	13
Sol-Vex®	14
Quick Reference Guide	15

Also available
Surgical Glove Catalog



1-800-952-9916
www.ansellhealthcare.com

Everything you touch...we touch.

At Ansell Healthcare, we understand, anticipate and respond to what you need in examination gloves. After all...gloves are your link and point of contact with your patients.

You asked for gloves that not only provide added protection, but also suit a variety of needs. We responded with a full line of premium examination gloves that enable you to perform an unlimited scope of procedures. Each glove is designed to meet your needs and all of our gloves offer comfort and reliable barrier protection. Whether you need powder-free latex, latex-free or something in between, we have the right examination gloves for you!

Ansell's Commitment to Innovation

Ansell examination gloves are manufactured with the highest level of quality available. And because we pledge to bring you, our customer, products with the latest innovations in the gloving market, you can always look to Ansell for the right product to fit your needs. Some of our latest innovations include:

Allergy Protection

Ansell's newest exam glove, Micro-Touch® NitraFree™, is manufactured without latex and sulfur-based chemical accelerators, helping to protect against Type I (latex) allergies and Type IV (chemical) allergies. Visit AnsellProtects.com for more information on allergy solutions.

Handfriendly Technology

Our handfriendly technology is available on the Micro-Touch® NextStep®. This latex glove offers a glycerol coating which hydrates, moisturizes and nourishes the skin.

Low-Protein Technology

Ansell utilizes advanced technology to reduce total water-based extractable protein thereby greatly reducing residual protein levels on latex gloves.

Chemotherapy Approved

Several of our exam gloves have been tested and are FDA approved for handling chemotherapy drugs.¹

Innovative Packaging

Ansell's newly designed exam glove pack is a reflection of our commitment to environmental sustainability. This innovative packaging features more exam gloves per box and per case in order to decrease packaging waste, increase storage space at your facility and reduce re-stocking orders to save on fuel consumption and shipping costs. Plus, the newly designed box will continue to fit into your existing wall brackets.



Look for our eco-friendly packaging logo to identify these products.



200 Pieces
Per Box

1. Specific chemical resistance data available upon request: 1-800-952-9916 (U.S.A. only). Gloves used for protection against chemotherapy drugs must be selected specifically for the type of chemicals used. Review material safety data sheets for the chemicals being used to determine the required level of protection.

Ansell's Value Proposition

We are committed to protecting you and your patients with the highest quality examination gloves. Here's how we do it.

Our Products

Our full line of premium examination gloves offers a wide variety of products to match the demanding requirements of the healthcare industry.

- Ansell has a rich heritage in product innovation, with over a century of industry milestones and manufacturing breakthroughs to our credit.
- Product concepts and performance are validated by internal laboratory testing. New product ideas are put through a rigorous Stage Gate® process to ensure that products function at expected performance levels. Healthcare professionals, in the appropriate healthcare environments, conduct final product evaluations and validations to ensure that our products stand up in real-life situations.
- Ansell is one of the first glove manufacturers to adopt a Six Sigma approach to operational excellence ensuring across-the-board efficiencies and quality.
- Ansell has implemented World Class Manufacturing (WCM) as a holistic initiative toward operational excellence.
- All of Ansell's examination gloves meet or exceed ASTM examination glove standards.

Our Service

We don't just sell products; we provide solutions to your specific needs.

- Our highly trained staff of sales representatives, clinical consultants and customer service department are without equal in the industry.
- Ansell addresses your needs with practical and actionable solutions.
- We offer a full range of educational programs designed to keep your staff up-to-date on the latest glove technology, trends and healthcare issues. Available offline or online, we are continuously providing new programs.
- Our one-of-a-kind 10-Step Glove Evaluation Process is the most thorough and comprehensive in the industry.
- We strive to always have the right product in your hands at the right time.

Our Value

Some customers require creative solutions to short cycle replenishment. Ansell partners with its network of distributors to ensure efficient distribution.

- Our Glove Management Program helps you streamline your inventory while eliminating misuse and redundancy.
- At Ansell, quality and value go hand in hand – we will work with you to maximize your facility's financial goals.
- We can negotiate flexible corporate agreements to suit the complexities of different provider organizations.
- Ansell teams up with you to provide customized logistical solutions, ensuring that your facility's individual requests are addressed timely and efficiently.
- We offer long-term Group Purchasing Organization (GPO) relationships. We understand your special purchasing needs and provide value-added solutions.

Latex-Free, Powder-Free Examination Gloves

Nitrile



SIZE	PRODUCT #
X-Small	6034510
Small	6034511
Medium	6034512
Large	6034513
X-Large	6034514

100 Gloves/Dispenser Box – 10 Boxes/Case

Micro-Touch® NitraFree™ Pink Nitrile, Powder-Free Examination Gloves

Micro-Touch® NitraFree™ pink nitrile examination gloves are manufactured without sulfur-based chemical accelerators*, helping to protect you from Type I (latex) and Type IV (chemical) allergies. Its premium nitrile material is soft with more stretch for greater comfort, fit and excellent protection in all risk settings. Proceeds from the purchase of every case go to Susan G. Komen for the Cure® in support of breast cancer research and outreach.†

- Material: 100% nitrile; no natural rubber latex or sulfur-based chemical accelerators*
- Color: Pink
- Textured fingertip finish
- Inverted, 9.5" beaded cuff
- FDA approved for handling chemotherapy drugs1**
- Ambidextrous
- Nonsterile
- Single use



*Does not contain any sulfur-based accelerators such as thiurams, carbamates, thiazoles, guanidine or thioureas.

†Through March 31, 2013, Ansell will donate \$2 per case sold of the Micro-Touch® NitraFree™ to Susan G. Komen for the Cure®, with a minimum guaranteed donation of \$25,000.

**Tested for use with the following 9 chemotherapy drugs: Vincristine Sulfate, Cyclophosphamide, Doxorubicin Hydrochloride, 5-Fluorouracil, Methotrexate, Etoposide, Paclitaxel, ThioTEPA and Carmustine.2

2. Caution: Testing showed an average breakthrough time of 45.3 minutes with Carmustine.



SIZE	PRODUCT #
X-Small	6034300
Small	6034301
Medium	6034302
Large	6034303
X-Large	6034304

200 Gloves/Dispenser Box – 10 Boxes/Case

Micro-Touch® Nitrile Nitrile, Powder-Free Examination Gloves

Micro-Touch® Nitrile examination gloves feature an advanced stretch formulation providing you with ultimate comfort and fit. With its high strength and excellent protection features, this glove is ideal for use in all risk settings. Now featured in a 200 count box to help you go GREEN. With more gloves in a box it means less packaging waste.

- Material: 100% nitrile
- Color: Blue
- Textured fingertip finish
- Inverted, 9.5" beaded cuff
- FDA approved for handling chemotherapy drugs1*
- Ambidextrous
- Nonsterile
- Single use




*Tested for use with the following 7 chemotherapy drugs: Dacarbazine, Cyclophosphamide, Doxorubicin Hydrochloride, 5-Fluorouracil, Cisplatin, Etoposide and Paclitaxel.

1. Specific chemical resistance data available upon request: 1-800-952-9916 (U.S.A. only). Gloves used for protection against chemotherapy drugs must be selected specifically for the type of chemicals used. Review material safety data sheets for the chemicals being used to determine the required level of protection.

Latex-Free, Powder-Free Examination Gloves

Nitrile



SIZE	PRODUCT #
X-Small	6034100
Small	6034101
Medium	6034102
Large	6034103
X-Large	6034104

100 Gloves/Dispenser Box – 10 Boxes/Case

Micro-Touch® Smooth Nitrile

Nitrile, Powder-Free Examination Gloves

Micro-Touch® Smooth Nitrile examination gloves provide the protection to meet the daily needs of the healthcare professional. Its anti-slip, smooth finish provides excellent tactile sensitivity, and its blue color clearly identifies it as a nitrile glove.

- Material: 100% nitrile
- Color: Blue
- Smooth finish
- Inverted, 9.5" beaded cuff
- FDA approved for handling chemotherapy drugs^{1*}
- Ambidextrous
- Nonsterile
- Single use



^{1*}Tested for use with the following 9 chemotherapy drugs: Carmustine, Cisplatin, Cyclophosphamide, Dacarbazine, ThioTEPA, 5-Fluorouracil, Doxorubicin Hydrochloride, Etoposide and Paclitaxel.



SIZE	PRODUCT #
X-Small	6034050
Small	6034051
Medium	6034052
Large	6034053
X-Large	6034054

100 Gloves/Dispenser Box – 10 Boxes/Case

Micro-Touch® Nitrile E.P.®

Nitrile, Powder-Free Examination Gloves

The Micro-Touch® Nitrile E.P.® examination gloves have an extended cuff that provides extra protection over the wrist and arm for healthcare professionals in high-risk clinical settings. The extra protection is beneficial when dealing with chemotherapy agents, bodily fluids and blood-borne pathogens. It provides excellent chemical and puncture resistance and is textured for better grip in wet and dry procedures.

- Material: 100% nitrile
- Color: Blue
- Textured fingertip finish
- Extended, 12" beaded cuff
- FDA approved for handling chemotherapy drugs^{1*}
- Ambidextrous
- Nonsterile
- Single use



^{1*}Tested for use with the following 9 chemotherapy drugs: Carmustine, Cisplatin, Cyclophosphamide, Dacarbazine, ThioTEPA, 5-Fluorouracil, Doxorubicin Hydrochloride, Etoposide and Paclitaxel.

1. Specific chemical resistance data available upon request: 1-800-952-9916 (U.S.A. only). Gloves used for protection against chemotherapy drugs must be selected specifically for the type of chemicals used. Review material safety data sheets for the chemicals being used to determine the required level of protection.

Latex-Free, Powder-Free Examination Gloves

Neoprene



SIZE	PRODUCT #
X-Small	3770
Small	3771
Medium	3772
Large	3773
X-Large	3774

100 Gloves/Dispenser Box – 10 Boxes/Case

*Tested for use with the following 17 chemotherapy drugs: Carmustine, Cisplatin, Doxorubicin Hydrochloride, Mitomycin-C, Melphalan, Carboplatin, Etoposide, Paclitaxel, 5-Fluorouracil, Daunorubicin Hydrochloride, Idarubicin Hydrochloride, Bleomycin Sulfate, ThioTEPA, Cyclophosphamide, Mitoxantrone, Ifosfamide and Cytarabine Hydrochloride.

Micro-Touch® Affinity®


Polychloroprene, Powder-Free Examination Gloves

Micro-Touch® Affinity® examination gloves provide the fit, comfort, tactile sensitivity and stretchability of natural rubber latex (NRL) without the risks associated with NRL or glove powder complications. A polymer coating provides excellent donnability for this polychloroprene glove.

- Material: 100% polychloroprene (neoprene)
- Color: Green
- Textured fingertip finish
- Inverted, 9.5" beaded cuff
- FDA approved for handling chemotherapy drugs^{1*}
- Ambidextrous
- Nonsterile
- Single use



Vinyl



SIZE	PRODUCT #
Small	3091
Medium	3092
Large	3093
X-Large	3094

100 Gloves/Dispenser Box – 10 Boxes/Case

Micro-Touch® Elite®


Polyvinyl Chloride, Powder-Free Examination Gloves

Micro-Touch® Elite® powder-free, stretch vinyl examination gloves are produced with a polyvinyl chloride formulation that provides a more comfortable fit than traditional vinyl gloves. They are an excellent alternative for latex-sensitive individuals for use in low-risk clinical applications that do not involve blood-borne pathogens or bodily fluids. This glove is not formulated with DEHP.

- Material: Stretch polyvinyl chloride (PVC)
- Color: Cream
- Smooth finish
- Inverted, 9.5" beaded cuff
- Ambidextrous
- Nonsterile
- Single use



Sterile, Powder-Free Examination Gloves



Micro-Touch® Nitra-Tex® Sterile Pairs and Singles

Sterile Nitrile, Powder-Free Examination Gloves

Micro-Touch® Nitra-Tex® Sterile exam gloves feature a 12" extended cuff and textured fingertips, providing superior barrier protection and grip in a sterile latex-free exam glove.

- Material: 100% nitrile
- Color: Blue
- Textured fingertip finish
- Extended, 12" beaded cuff
- FDA approved for handling chemotherapy drugs^{1*}
- Ambidextrous
- Sterile (for nonsurgical procedures)
- Single use



SIZE	PAIRS PRODUCT #	SINGLES PRODUCT #
Small	6034151	6034251
Medium	6034152	6034252
Large	6034153	6034253
X-Large	6034154	6034254

PAIRS:
50 Pairs/Dispenser Box – 4 Boxes/Case

SINGLES:
100 Gloves/Dispenser Box – 4 Boxes/Case

*Tested for use with the following 8 chemotherapy drugs: Vincristine Sulfate, Cyclophosphamide, Doxorubicin Hydrochloride, 5-Fluorouracil, Methotrexate, Etoposide, Paclitaxel and ThioTEPA.

1. Specific chemical resistance data available upon request: 1-800-952-9916 (U.S.A. only). Gloves used for protection against chemotherapy drugs must be selected specifically for the type of chemicals used. Review material safety data sheets for the chemicals being used to determine the required level of protection.



Micro-Touch® Plus Sterile Singles

Sterile Latex, Powder-Free Examination Gloves

Micro-Touch® Plus Sterile powder-free, textured, natural rubber latex examination gloves feature a polymer coating to aid in donning and provide excellent tactile sensitivity, durability and grip that healthcare professionals can expect in both wet and dry procedures. This sterilized natural rubber latex exam glove is recommended for sterile exam (non-surgical) procedures. Advanced manufacturing processes reduce total water extractable protein to 50 µg or less per gram.

- Material: Natural rubber latex
- Color: Cream
- Fully textured finish
- Inverted, 9.5" beaded cuff
- Ambidextrous
- Sterile (for nonsurgical procedures)
- Single use



SIZE	PRODUCT #
Small	6016001
Medium	6016002
Large	6016003

100 Gloves/Dispenser Box – 4 Boxes/Case

Caution: Micro-Touch® Plus Sterile Singles contain natural rubber latex, which may cause allergic reactions. Safe use of this glove by or on latex-sensitized individuals has not been established.

Latex, Powder-Free Examination Gloves



SIZE	PRODUCT #
X-Small	6015300
Small	6015301
Medium	6015302
Large	6015303
X-Large	6015304

150 Gloves/Dispenser Box – 10 Boxes/Case

Micro-Touch® Plus


Latex, Powder-Free Examination Gloves

The Micro-Touch® Plus examination glove is made from natural rubber latex and features a polymer coating for ease of donning, with a textured wet-grip surface. It is powder-free to help eliminate powder-related complications and is manufactured to meet or exceed ASTM standards. Advanced manufacturing processes reduce total water-extractable protein to 50 µg or less per gram. Now featured in a 150 count box to help you go **GREEN**. With more gloves in a box it means less packaging waste.

- Material: Natural rubber latex
- Color: Cream
- Fully textured finish
- Inverted, 9.5" beaded cuff
- Ambidextrous
- Nonsterile
- Single use



Eco-friendly packaging



SIZE	PRODUCT #
X-Small	3200
Small	3201
Medium	3202
Large	3203
X-Large*	3204

100 Gloves/Dispenser Box – 10 Boxes/Case

*80 Gloves/Dispenser Box – 10 Boxes/Case

Micro-Touch® NextStep®

Latex, Powder-Free, Handfriendly Examination Gloves

Micro-Touch® NextStep® gloves with a glycerol moisturizer provide superior protection and tactile sensitivity, along with the comfort and advantages of latex. Designed with the healthcare professional's daily work routine in mind, the light green exam gloves are easy to don, are textured for reliable grip, and nourish skin during use. Advanced manufacturing processes reduce total water-extractable protein to 50 µg or less per gram.

- Material: Natural rubber latex
- Color: Light green
- Fully textured finish
- Inverted, 9.5" beaded cuff
- Ambidextrous
- Nonsterile
- Single use



Latex, Powder-Free Examination Gloves



SIZE	PRODUCT #
X-Small	9201
Small	9203
Medium	9205
Large	9207

1 Pair/Package – 25 Pairs/Box – 4 Boxes/Case

Micro-Touch® E.P.® Gloves

Latex, Powder-Free Examination Gloves

Micro-Touch® E.P.® Gloves are made from natural rubber latex, feature extra long reinforced cuffs, and are thicker for extra protection, making them ideal for use in higher-risk clinical applications. These powder-free gloves are hand specific which allow for added comfort and fit.

- Material: Natural rubber latex
- Color: Cream
- Fully textured finish
- Extended straight, 12" reinforced cuff
- FDA approved for handling chemotherapy drugs^{1*}
- Hand specific
- Nonsterile
- Single use



* Tested for use with the following 18 chemotherapy drugs: Bleomycin Sulfate, Carboplatin, Carmustine, Cisplatin, Cyclophosphamide, Cytarabine Hydrochloride, Daunorubicin Hydrochloride, Doxorubicin, Etoposide, 5-Fluorouracil, Idarubicin, Ifosfamide, Mechlorethamine Hydrochloride, Melphalan, Mitomycin, Mitoxantrone, Paclitaxel and ThioTEPA.

Latex, Powdered Examination Gloves



SIZE	PRODUCT #
Small	6116301
Medium	6116302
Large	6116303

150 Gloves/Dispenser Box – 10 Boxes/Case

Micro-Touch® Powder Latex, Powdered Examination Gloves

The Micro-Touch® Powder glove features a smooth finish for better tactile sensitivity and a straight reinforced cuff. This glove is made from natural rubber latex for superior hand barrier protection and is powdered with absorbable dusting powder, U.S.P. These gloves are manufactured utilizing Ansell's unique ergonomic and concave formers, which provide unparalleled fit. Now featured in a 150 count box to help you go **GREEN**. With more gloves in a box it means less packaging waste.

- Material: Natural rubber latex
- Color: Cream
- Smooth finish
- Straight, 9.5" reinforced cuff
- Ambidextrous
- Nonsterile
- Single use



Eco-friendly packaging



General Purpose Gloves



SIZE	PRODUCT #
Small	8984
Medium	8986
Large	8988
X-Large	8980

12 Pairs/Bag – 12 Bags/Case

Housekeeping Latex, Powder-Free Gloves

For gloves that stand up to the rigors of everyday routine cleaning and sanitizing, to handling chemicals and harsh materials, Ansell offers alternatives to keep hands protected. Our powder-free housekeeping glove features a soft cotton flock lining for comfort and easy donnability. This glove provides an embossed, sure-grip surface and dependable strength to ensure that your staff is protected.

- Material: Natural rubber latex
- Color: Yellow
- Embossed, sure-grip surface
- 12" length
- Reusable



SIZE	PRODUCT #
Small	5789015
Medium	5789016
Large	5789017
X-Large	5789018

12 Pairs/Bag – 4 Bags/Case

Latex/Nitrile Utility Latex, Powder-Free Gloves

For gloves that stand up to the rigors of everyday routine cleaning and sanitizing, to handling chemicals and harsh materials, Ansell offers alternatives to keep hands protected. Our powder-free utility glove protects against acids, oils, disinfectants, detergents and abrasives. They also provide a diamond-embossed, sure-grip surface and dependable strength to ensure that your staff is protected.

- Material: Natural rubber latex and nitrile blend
- Color: Blue
- Diamond-embossed, sure-grip surface
- 13" length
- Reusable



General Purpose Gloves



SIZE	PRODUCT #
6-6.5	117272
7-7.5	117273
8-8.5	117274
9-9.5	117275
10	117276
11	117277

12 Pairs/Bag – 12 Bags/Case

Sol-Vex®

Nitrile, Powder-Free Gloves

Sol-Vex® Nitrile, powder free gloves are crafted from a versatile, high-functioning polymer to create an ideal balance between protection and performance. Regardless of the environment (wet or dry) these gloves will provide both strength and chemical protection, without compromising on flexibility, comfort or dexterity.

- Material: 100% nitrile
- Color: Green
- Sandpatch embossed, sure grip surface
- Chemical resistant
- 13" length
- Reusable



Product	Size	Product #	Page #	Material	Powder	Color	Beaded Cuff	Former	Finish	Minimum Length (in)	Fingertip – Single Thickness (mils)	Case Configuration
Latex-Free, Powder-Free Examination Gloves												
Nitrile												
Micro-Touch® NitraFree™	XS	6034510	6	Nitrile	No	Pink	Yes	Straight	Textured Fingertips	9.4	3.1	100 Gloves/Box 10 Boxes/Case
	S	6034511										
	M	6034512										
	L	6034513										
	XL	6034514										
Micro-Touch® Nitrile	XS	6034300	6	Nitrile	No	Blue	Yes	Straight	Textured Fingertips	9.4	3.1	200 Gloves/Box 10 Boxes/Case
	S	6034301										
	M	6034302										
	L	6034303										
	XL	6034304										
Micro-Touch® Smooth Nitrile	XS	6034100	7	Nitrile	No	Blue	Yes	Straight	Smooth	9.4	3.9	100 Gloves/Box 10 Boxes/Case
	S	6034101										
	M	6034102										
	L	6034103										
	XL	6034104										
Micro-Touch® Nitrile E.P.®	XS	6034050	7	Nitrile	No	Blue	Yes	Straight	Textured Fingertips	11.6	4.7	100 Gloves/Box 10 Boxes/Case
	S	6034051										
	M	6034052										
	L	6034053										
	XL	6034054										
Neoprene												
Micro-Touch® Affinity®	XS	3770	8	Polychloroprene (neoprene)	No	Green	Yes	Straight	Textured Fingertips	9.4	6.3	100 Gloves/Box 10 Boxes/Case
	S	3771										
	M	3772										
	L	3773										
	XL	3774										
Vinyl												
Micro-Touch® Elite®	S	3091	8	Stretch Polyvinyl Chloride	No	Cream	Yes	Straight	Smooth	9.1	3.9	100 Gloves/Box 10 Boxes/Case
	M	3092										
	L	3093										
	XL	3094										
Sterile, Powder-Free Examination Gloves												
Micro-Touch® Nitra-Tex® Sterile Pairs	S	6034151	9	Nitrile	No	Blue	Yes	Straight	Textured Fingertips	11.6	4.7	50 Pairs/Box 4 Boxes/Case
	M	6034152										
	L	6034153										
	XL	6034154										
Micro-Touch® Nitra-Tex® Sterile Singles	S	6034251	9	Nitrile	No	Blue	Yes	Straight	Textured Fingertips	11.6	4.7	100 Gloves/Box 4 Boxes/Case
	M	6034252										
	L	6034253										
	XL	6034254										
Micro-Touch® Plus Sterile Singles	S	6016001	9	Natural Rubber Latex	No	Cream	Yes	Straight	Textured	9.4	5.5	100 Gloves/Box 4 Boxes/Case
	M	6016002										
	L	6016003										
Latex, Powder-Free Examination Gloves												
Micro-Touch® Plus	XS	6015300	10	Natural Rubber Latex	No	Cream	Yes	Straight	Textured	9.1	3.1	150 Gloves/Box 10 Boxes/Case
	S	6015301										
	M	6015302										
	L	6015303										
	XL	6015304										
Micro-Touch® NextStep®	XS	3200	10	Natural Rubber Latex	No	Light Green	Yes	Straight	Textured	9.5	7.1	100 Gloves/Box 10 Boxes/Case *XL – 80 Gloves/Box 10 Boxes/Case
	S	3201										
	M	3202										
	L	3203										
	XL*	3204										
Micro-Touch® E.P.® Gloves	XS	9201	11	Natural Rubber Latex	No	Cream	No	Curved	Textured	12.0	9.1	1 Pair/Package 25 Pairs/Box 4 Boxes/Case
	S	9203										
	M	9205										
	L	9207										
Latex, Powdered Examination Gloves												
Micro-Touch® Powder	S	6116301	12	Natural Rubber Latex	Yes	Cream	No	Straight	Smooth	9.1	6.7	150 Gloves/Box 10 Boxes/Case
	M	6116302										
	L	6116303										
General Purpose Gloves												
Housekeeping	S	8984	13	Natural Rubber Latex	No	Yellow	No	Curved	Embossed	12.0	22	12 Pairs/Bag 12 Bags/Case
	M	8986										
	L	8988										
	XL	8980										
Latex/Nitrile Utility	S	5789015	13	Natural Rubber Latex/ Nitrile	No	Blue	No	Curved	Embossed	13.0	20	12 Pairs/Bag 4 Bags/Case
	M	5789016										
	L	5789017										
	XL	5789018										
Sol-Vex®	6-6.5	117272	14	Nitrile	No	Green	No	Curved	Embossed	13.0	15	12 Pairs/Bag 12 Bags/Case
	7-7.5	117273										
	8-8.5	117274										
	9-9.5	117275										
	10	117276										
	11	117277										

**For additional information and samples, call Ansell
at 800-952-9916 or visit us at www.ansellhealthcare.com**

Ansell Healthcare Products LLC | ansellhealthcare.com
200 Schulz Drive | Red Bank, NJ 07701 U.S.A. | Phone: (800) 952-9916 (U.S.A. only)

Ansell, Micro-Touch, NitraFree, E.P., Elite, Affinity, NitraTex and NextStep are trademarks owned by Ansell Limited or one of its affiliates.
©2012 Ansell Healthcare Products LLC. All Rights Reserved.
The COLOR PINK® and PINK NITRILE® are registered trademarks owned by Colur World, LLC. 2006 All Rights Reserved.

BR140 8/2012

Printed on recycled paper using soy inks.



Ansell
Everything you touch...we touch.