

NEW!

PowerSpeed™

STRONGER THAN POWERFLEX

THE ULTIMATE SPEED TAPE

- Designed/developed for speed taping or for a stronger spat
- 30% stronger than standard PowerFlex
- Delivers maximum support in a modified taping system
- Won't slip from sweat or water
- Hand-tear



Speed Taping (modified taping) to the Ankle:

- Apply PowerSpeed as your base layer
- Apply figure eights and heel locks, 3 short stirrups
- Finish with figure eights and heel lock



STEP 1:

Apply PowerSpeed directly to the skin - removing all stretch - for 32 lbs. of tensile strength.



STEP 2:

Apply athletic tape or PowerTape over PowerSpeed for up to 72 lbs. of tensile strength.

ORDERING INFORMATION

U.S. Patent No. 5,762,623 European Patent No. 0839222

Catalog#	Size	Rolls/Case
PS3720XX-024	2" x 6 yds.	24
PS3725XX-016	2 3/4" x 6 yds.	16
PS3730XX-016	3" x 6 yds.	16

Colors Available: Black (BK), White (WH).

When placing an order, please state catalog # and color. Example: PS3720BK. XX = Color Code.



Andover Healthcare, Inc.

9 Fanaras Drive, Salisbury, MA 01952 Toll Free: 800.432.6686 Phone: 978.465.0044 Fax: 978.462.0003

www.andoverhealthcare.com

Proudly Made in USA