



◀ 1 of 3 ▶

[Product Advantages](#)

Piston Syringes and Trays

Piston Syringes

70cc piston syringe features double-flanged plunger tip to ensure adequate vacuum, as well as a thumb ring and smooth piston action for easy one-hand use. Single-use, sterile.

Piston (Toomey)

70cc piston syringe barrel fits all standard cystoscopes or resectoscopes. Single-use, sterile

Order No.	Description	Pkg
0038470	70cc piston irrigation syringe	50/case
750301	Piston syringe, luer adapter and graduated flask packaged in basin with CSR overwrap	20/case
750307	Piston syringe, luer adapter, graduated flask and underpad packaged in basin with peelback lid	20/case
750310	Piston syringe and graduated flask packaged in basin with peelback lid	20/case
0038460	Toomey syringe with catheter and luer adapter tip	50/case