

BD PhaSeal™ System

IV Bag Access Device

Work safe. Enjoy life.



Infusion Adapter (C100)



The Infusion Adapter has a built-in BD PhaSeal™ Connector to facilitate closed addition to IV bags. The adapter allows the user the opportunity to select the tubing they need to spike into the bag. A non-vented IV set should be used together with the Infusion Adapter. If a vented set is used, make sure that the air inlet is closed throughout the infusion procedure. PVC- and BPA-free.



Helping all people
live healthy lives

Physical characteristics

REF / ART. NO.	515306
LENGTH	90 mm (3.54 in.)
WIDTH	24 mm (0.94 in.)
HEIGHT	43 mm (1.69 in.)
INTERNAL VOLUME	0.16 mL
PACKAGING: USER BOX	50 pcs
PACKAGING: TRANSPORT BOX	200 pcs