

# BD PhaSeal™ System

## Syringe Safety Device

Work safe. Enjoy life.



### Injector Luer Lock (N35)



The Injector Luer Lock ensures a closed drug transfer by means of double, elastomeric membranes. The Injector Luer Lock attaches to a standard syringe or IV tubing for a dry, leakproof connection during preparation and administration.



Helping all people  
live healthy lives

#### Physical characteristics

REF / ART. NO.	515003
EXTERNAL FITTING	Luer Lock
DIMENSION LENGTH	71 mm (2.79 in.)
DIMENSION WIDTH	17 mm (0.67 in.)
CANNULA VOLUME	0.04 mL
PACKAGING: USER BOX	50 pcs
PACKAGING: TRANSPORT BOX	200 pcs