

BD PhaSeal™ System

Drug Vial Access Devices

Work safe. Enjoy life.



Protector



The Protector is a pressure equalization device, permanently attached to the vial used in the preparation of drugs from either powder or liquid form. When attached to a vial, the Protector provides closed pressure equalization through the sealed expansion chamber and dry connections via the double membrane design. This ensures that there is neither overpressure nor vacuum when air or fluid is injected into or aspirated from the vial.



Helping all people
live healthy lives

Physical characteristics

PRODUCT NAME	Protector (P14)	Protector (P21)	Protector (P28)	Protector (P50)
REF / ART. NO.	515100	515102	515104	515105
FITS TO VIAL NECK SIZE	013 mm	020 mm	028 mm	020 mm
PRESSURE EQUALIZATION	20 mL	20 mL	60 mL	50 mL
COLOR PROTECTIVE CAP	Blue ■	Red ■	White □	Green ■
PACKAGING: USER BOX	50 pcs	50 pcs	25 pcs	25 pcs
PACKAGING: TRANSPORT BOX	200 pcs	200 pcs	100 pcs	100 pcs

BD Medical
1 Becton Drive
Franklin Lakes, NJ 07417
www.bd.com/phaseal