

Venous Blood Collection

BD Vacutainer® Passive Shielding
Blood Collection Needle



BD


Helping all people
live healthy lives


*A New Generation of Safety from BD,
the leader in safety-engineered technology*

The BD Vacutainer® Passive Shielding Blood Collection Needle is designed to provide the highest amount of protection against needlesticks during a blood draw. The needle's safety shield releases automatically upon insertion of the first tube, providing protection from the time the first tube is inserted through product disposal.

BD Vacutainer® Passive Shielding Blood Collection Needle

Ordering Information

368636	Passive Shielding Blood Collection Needle	21g x 1" needle	25/Box 200/Case	
--------	--	--------------------	--------------------	---

368637	Passive Shielding Blood Collection Needle	22g x 1" needle	25/Box 200/Case	
--------	--	--------------------	--------------------	---



A part of the BD Vacutainer® Family of Venous Blood Collection Products

CAUTION: Handle all biologic samples and blood collection "sharps" (lancets, needles, luer adapters, and blood collection sets) in accordance with the policies and procedures of your facility. Obtain appropriate medical attention in the event of any exposure to biologic samples (e.g., through a puncture injury) since samples may transmit viral hepatitis, HIV (AIDS), or other infectious diseases. Utilize any safety-engineered feature if the blood collection device provides one. Discard all blood collection "sharps" in biohazard containers approved for their disposal.

BD Global Technical Services:

1.800.631.0174

BD Customer Service:

1.888.237.2762

www.bd.com/vacutainer



BD Diagnostics

Prenalytical Systems

1 Becton Drive

Franklin Lakes, NJ 07417

BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company.

©2005 BD. Printed in USA 01/05 VS7355