

ICON MONO

ONE-STEP ANTIBODY TEST FOR INFECTIOUS
MONONUCLEOSIS



When all the symptoms point to 'IM', use ICON Mono to be sure.

ICON Mono from Beckman Coulter provides fast, accurate and easy-to-read results at the point-of-care.

Blood Banking
Centrifugation
Chemistry
Flow Cytometry
Hematology
Immunoassay
Information Systems
Lab Automation
Molecular Diagnostics
Rapid Diagnostics

When your patients present with fever, sore throat, swollen lymph glands or other symptoms common with a viral infection, you may not be able to distinguish these from infectious mononucleosis (IM).

That's why there's ICON Mono.

- Combines finger-stick whole blood convenience with EIA membrane technology for fast, accurate results.
- Uses direct solid-phase immunoassay technology for early qualitative detection of IM heterophile antibodies in whole blood, human serum or plasma.
- 99.9% relative sensitivity and 98.8% relative specificity serum/plasma; 98.3% relative specificity WB assay.

Make it a point to use the one-step, easy-to-read rapid test proven to deliver fast and accurate results – ICON Mono from Beckman Coulter.

- Highly accurate – 99.9% relative sensitivity and 98.8% relative specificity serum/plasma; 98.3% relative specificity WB assay.
- Single reagent procedure provides easy-to-read results in 3-8 minutes.
- Walk-away procedure accommodates interruptions in workflow.
- Built-in procedural quality control ensures proper test performance.
- Whole-blood procedure allows for specimen to be introduced directly into the test well – no messy pretreatment in separate tube.
- CLIA-waived status for whole blood.

CLIA

- CPT Code* #86308QW, 86308
- Whole Blood - Waived
- Serum/Plasma - Moderately Complex



Ordering Information

ICON Mono Test Kit

395100

- 30 Test Cassettes
- Developer Solution
- Negative Control
- Positive Control
- 35 Capillary Tubes
- Bulb

Beckman Coulter is the leader in fecal occult blood testing and offers the very best in rapid tests designed to simplify the diagnostic process. Early detection is the single most important advantage a patient can have to survive life-threatening diseases. Beckman Coulter provides the tools and tests physicians need to aid in the diagnosis and management of a variety of disease states.

Beckman Coulter also provides a variety of Immunoassay, Chemistry, Hematology and Centrifugation instruments to aid in the diagnosis and management of a full spectrum of medical conditions. Each system's compact design and cost-effective performance makes it a perfect complement to a physician's office laboratory.

To order product, contact your medical supply distributor. For technical assistance regarding rapid test kits, call Technical Marketing at 800-877-6242 or email askpcd@beckman.com.

*Exemplar CPT provided. Refer to regional CMS fee schedule for relevant/current CPT codes and reimbursement rates.

Lab Testing Made Easy

Beckman Coulter Rapid Diagnostics Product Offering

Hemoccult

Hemoccult SENSE

Hemoccult II SENSE elite

Hemoccult ICT

CLIA-waived rapid tests for the detection of hidden blood in feces

Gastroculta

CLIA-waived rapid test for detecting hidden blood in gastric fluid and determining gastric pH

ICON Hp

CLIA-waived (whole blood) and moderately complex (serum) rapid test for the detection of antibodies to *H. pylori*

ICON 20 hCG

ICON 25 hCG

CLIA-waived (urine) and moderately complex (serum) rapid tests for the detection of hCG

ICON Mono

CLIA-waived (whole blood) and moderately complex (serum) antibody test for infectious mononucleosis

ICON DS Strep A

ICON SC Strep A

CLIA-waived dipstick and cassette rapid format tests for the detection of Strep A, using throat swab sample



Beckman Coulter, the stylized logo, Hemoccult, Hemoccult II, SENSE, Gastroculta and ICON are trademarks of Beckman Coulter, Inc. and are registered with the USPTO.

For Beckman Coulter's worldwide office locations and phone numbers, please visit www.beckmancoulter.com/contact

DS-17652A B2012-13776

www.beckmancoulter.com

© 2013 Beckman Coulter, Inc.