

# ICON<sup>SC</sup> Strep A

DIRECT GROUP A STREPTOCOCCAL ANTIGEN TEST



**Get your Strep A patients on the move.**

**Get fast, accurate results with this simple, easy test.**



ICON SC Strep A is a simple, easy-to-use, CLIA-waived, cassette-format test. In as little as 5 minutes, you will have accurate results so you can begin your patient's treatment immediately and confidently. With outstanding performance, even in low-colony-count samples and simple-to-interpret test results, there is little chance for misdiagnosis, minimizing the chance of rheumatic fever or overprescription of antibiotics.

- **98% accuracy**
- **Results in 5 minutes**
- **Easy-to-use cassette format**
- **On-board extraction**
- **Fast, easy and accurate)**
- **Storage at 2-30°C (35-86°F)**

Blood Banking  
Centrifugation  
Chemistry  
Flow Cytometry  
Hematology  
Immunoassay  
Information Systems  
Lab Automation  
Molecular Diagnostics  
**Rapid Diagnostics**





## PERFORMANCE DATA: MIXED POPULATION

In clinical trials, test sensitivity among the pediatric population was 96.2% (127/132).  
Test specificity was 98.7% (368/373) and accuracy was 98% (495/505).

## COMPETITIVE COMPARISON

	2	1	3
	Old ICON SC	NEW ICON SC	Quidel Quick View
Sensitivity	90.5%	96.2%	92%
Specificity	94.5	98.7	99
Accuracy	93.1	98	96

Data based on Beckman Coulter, Inc. and Quidel Corporation controlled tests

## CLIA

- CPT Code\* #87880QW - Strep A Test
- Waived



## Ordering Information

### ICON SC Strep A Test Kit

395161

- 25 Test Cassettes
- 25 Throat Swabs
- 1 Vial Reagent A
- 1 Vial Reagent B
- 1 Bottle Positive Control
- 1 Bottle Negative Control

\*Exemplar CPT provided. Refer to regional CMS fee schedule for relevant/current CPT codes and reimbursement rates.

Beckman Coulter is the leader in fecal occult blood testing and offers the very best in rapid tests designed to simplify the diagnostic process. Early detection is the single most important advantage a patient can have to survive life-threatening diseases. Beckman Coulter provides the tools and tests physicians need to aid in the diagnosis and management of a variety of disease states.

Beckman Coulter also provides a variety of Immunoassay, Chemistry, Hematology and Centrifugation instruments to aid in the diagnosis and management of a full spectrum of medical conditions. Each system's compact design and cost-effective performance makes it a perfect complement to a physician's office laboratory.

To order product, contact your medical supply distributor. For technical assistance regarding rapid test kits, call Technical Marketing at 800-877-6242 or email [askpcd@beckman.com](mailto:askpcd@beckman.com).

## Lab Testing Made Easy

### Beckman Coulter Rapid Diagnostics Product Offering

#### Hemoccult

#### Hemoccult SENA

#### Hemoccult II SENA elite

#### Hemoccult ICT

CLIA-waived rapid tests for the detection of hidden blood in feces

#### Gastroculta

CLIA-waived rapid test for detecting hidden blood in gastric fluid and determining gastric pH

#### ICON Hp

CLIA-waived (whole blood) and moderately complex (serum) rapid test for the detection of antibodies to *H. pylori*

#### ICON 20 hCG

#### ICON 25 hCG

CLIA-waived (urine) and moderately complex (serum) rapid tests for the detection of hCG

#### ICON Mono

CLIA-waived (whole blood) and moderately complex (serum) antibody test for infectious mononucleosis

#### ICON DS Strep A

#### ICON SC Strep A

CLIA-waived dipstick and cassette rapid format tests for the detection of Strep A, using throat swab sample



Beckman Coulter, the stylized logo, Hemoccult, Hemoccult II, SENA, Gastroculta and ICON are trademarks of Beckman Coulter, Inc. and are registered with the USPTO.

For Beckman Coulter's worldwide office locations and phone numbers, please visit [www.beckmancoulter.com/contact](http://www.beckmancoulter.com/contact)

DS-12540B B2012-13718

[www.beckmancoulter.com](http://www.beckmancoulter.com)

© 2013 Beckman Coulter, Inc.