



Product Description:

Bed Pan/Urinal Dispenser

- Holds a bed pan and up to two urinal bottles
- Keyholes for wall mounting
- Made of quartz powder coated aluminum
- Made in the USA

Dispenses:

- Most brand name standard bed pans and urinals

Product Specifications (Overall External Dimensions):

- 12 ⁵/₁₆" W x 12" H x 4 ¹¹/₁₆" D (31.3 cm x 30.5 cm x 11.8 cm)
- 1.9 lbs (0.8 kg) - approximated

Package Specifications:

- Package Qty 2 per case
- 16" W x 12" H x 12" D (40.6 cm x 30.5 cm x 30.5 cm)
- 6 lbs (2.7 kg) - approximated

Mounting Specifications:

