

AARON® 940

Bovie®



ELECTROSURGERY FOR THE OFFICE AND SURGI-CENTER



HIGH-FREQUENCY DESICCATOR

POWER

40 watt high frequency desiccator

FLEXIBILITY

Adjustable in 1/10th watt increments in the first 10 watts, in monopolar and bipolar

ACCURACY

More accurate power with Universal Power Control 100-240 VAC

SIMPLICITY

Less to do and less to buy makes this the most cost effective and user friendly unit on the market

COST EFFECTIVE

The Aaron® 940™ is a ground reference unit which means no grounding pad is required.

RELIABILITY

The Aaron® 940™ is backed by Bovie's four (4) year manufacturer's limited warranty.

AARON® 940™
INCLUDES

| | | |
|---|--|---|
| 3-Button Autoclavable Handpiece | Blunt & Sharp Dermal Tips | Wall Mount Kit |
|  |  |  |
| A901 | A806 & A804 (50 each) | A837 |

OTHER RECOMMENDED
ACCESSORIES

| | | |
|---|--|---|
| Footswitch | Reusable Bipolar Forceps* | Forcep Cable |
|  |  |  |
| A803 | A824 & A821 | A827BP |

| Disposable and Reusable Electrodes* | | | |
|---|---|---|---|
| Ball | Needle | Sterile Sheath | Non-Sterile Sheath |
|  |  |  |  |
| ES20 | ES02 | A910ST | A910 |


*More styles available

Industry-standard accessories fit most electrosurgical generators. Additional accessories available —contact your Bovie representative for details.

TECHNICAL
SPECIFICATIONS

Technical Specifications...

| | | |
|---|--|---------------------------|
| Dimensions: | Line Voltage: | Line Frequency: |
| Height 16.0 cm (6.3 in.) Width 22.8 cm (9 in.) Depth 11.4 cm (4.5 in.) Weight < 1.4 kg (< 3 lbs) | Universal Power Supply A940: 100 to 240 VAC | 50 - 60 Hz |
| Output Power: | Output Frequency: | Warranty: |
| To 40 watts (0-10 watts in 1/10th watt increments on both mono and bipolar) | 550 kHz | four (4) year (unit only) |



Bovie Medical Corporation • 5115 Ulmerton Road • Clearwater, FL 33760 • www.boviemed.com • sales@boviemed.com
U.S. Phone 1-800-537-2790 • Fax 1-800-323-1640 • Int'l Phone +1-727-384-2323 • Fax +1-727-347-9144

Bovie® is a registered trademark of Bovie Medical Corporation