


AARON® 950-G

*Bovie®*



# ELECTROSURGERY FOR THE OFFICE AND SURGERY CENTER



## COMPLETE

The Aaron 950-G is packaged with the essential accessories needed to perform an OB/GYN LEEP procedure.

## CONVENIENCE

The full system is situated on a convenient mobile stand for effortless portability.

## AFFORDABLE

The Aaron 950-G is an extremely affordable smoke evacuation system.

## PERFECTION


A variety of LLETZ loop electrodes are available for perfect surgical selection.

## RELIABILITY

The Aaron 950 unit is backed by Bovie's four (4) year manufacturer's limited warranty.

ELECTROSURGICAL GENERATOR

# AARON<sup>®</sup> 950™ -G INCLUDES

3-Button Autoclavable Handpiece	Non-Sterile Sheath	Blunt & Sharp Dermal Tips
		
A901	A910 (10 Included)	A806 & A804 (50 each)

Disposable Electrodes	Solid Disposable Grounding Pad	Reusable Grounding Cord	Wall Mount Kit
			
ES01 (3 Included), ES02, ES20	A1202 (5 Included)	A1202C	A837


Smoke Shark™ II Smoke Evacuator with 35-Hour Filter and Pneumatic Footswitch	6' Non-sterile Tube & Reducer Fitting	Mobile Stand & Tray
		
SE02	786T & SERF	A812-C

ESLK - Disposable Electrosurgical Loop Electrode Kit	ES07- Extended Ball Disposable Electrode - 5/box
	

Additional accessories available —contact your Bovie representative for details.

# TECHNICAL SPECIFICATIONS

## A950 Technical Specifications...

Dimensions:		Line Voltage:	Line Frequency:
Height	16.0 cm (6.3 in.)	A950: 120 VAC	50 - 60 Hz
Width	22.8 cm (9 in.)	A950-220: 240 VAC	
Depth	11.4 cm (4.5 in.)		
Weight	< 3.2 kg (< 7 lbs)		
Monopolar	Output Power	Output Frequency	Rep Rate
Cut	60 W @ 500 Ω	357 kHz ± 50 kHz	N/A
Blend	60 W @ 800 Ω	357 kHz ± 50 kHz	30 kHz ± 5 kHz
Coagulation	60 W @ 1000 Ω	375 kHz ± 50 kHz	60 kHz ± 5 kHz
Fulguration	35 W @ 1000 Ω	575 kHz ± 50 kHz	30 kHz ± 5 kHz
Bipolar	30 W @ 200 Ω	520 kHz ± 14 kHz + 50 kHz	19 kHz ± 5 kHz
<b>Warranty:</b> four (4) year on unit only 			

## SE02 Technical Specifications...

Dimensions: (without filter)	
Height	15.2 cm (6 in.)
Width	27.9 cm (11 in.)
Depth	39.4 cm (15.5 in.)
Weight	4.3 kgs (10 lbs)
Line Voltage:	
SE02: 120 VAC	
SE02-220: 240 VAC	
Line Frequency:	
50 - 60 Hz	
Noise Level, dBA: 55.0 dBA (Maximum)	
Warranty: limited two (2) year on unit only	