

# ProBarrier™ disposa chain™

## Disposable Bib Chain



***Eliminate germs, make-up, sweat and splatter with disposa chain!***

Designed to replace metal bib chains and clips. Use a fresh ProBarrier disposa chain for each patient and eliminate the need to reprocess bib chain.

### **Disposable:**

- Helps reduce the risk of cross-contamination. Use once and discard with each patient. Eliminates the need to disinfect metal bib chains.

### **Economical:**

- (300) 3/4" x 16-1/2" disposa chains per box. 15% more strips than the competition.

### **Patient Comfort:**

- Soft, lightweight non-latex material fits comfortably on patient. Eliminates cold chains that snag hair. Adhesive strips securely hold disposa chain to bib.



<b>ProBarrier disposa chain Description</b>	<b>Case Qty</b>	<b>Order#</b>
Disposable Bib Chain - 3/4" x 16" - 300/box	10	DC-0115S

## ProBarrier™ dispo chain™

### Disposable Bib Chain

#### Description:

ProBarrier dispo chain is a disposable bib chain.

#### Intended Use:

ProBarrier dispo chain is designed to replace metal bib chains and clips and to help prevent cross-contamination.

#### Product Specifications:

Material:	Non-latex material that will comfortably fit on patient and provide greater patient and user visibility.
Construction:	Adhesive strips on both ends.
Packaging:	Package is compact and fits into an operatory drawer.
Dimensions:	3/4" x 16-1/2".

#### Product Capabilities:

ProBarrier dispo chain saves money, improves patient comfort and reduces cross-contamination. ProBarrier dispo chain eliminates the labor and material costs associated with sterilizing as well as replacing metal chains and clips that are accidentally thrown away. ProBarrier dispo chain will not "chill" the patient's neck or snag hair. A fresh dispo chain used for each patient and discarded minimizes the chance for cross-contamination.

#### Accessory Products:

Other products available from Certol include:



**ProBarrier™ Chair Sox™**  
Environmental Surface  
Barriers



**ProBarrier™ Tube Sox™**  
Environmental Surface  
Barriers



**ProBarrier™ X-Ray Sox™**  
Environmental Surface  
Barriers



**ProTray™**  
Disposable Instrument  
Trays