

Optimizer[®] Polyp Trap



Multiple Chamber Polyp Trap

Applied Ingenuity[™] |

 **CONMED[™]**
ENDOSCOPIC TECHNOLOGIES

Optimizer[®] Polyp Trap

Optimizer[®] Key Features

VERSATILE

- Four independent, selectable trap chambers
- Multiple direct bypass channels
- Visual markers for each chamber

FAST


- Convenient central vacuum lumen

INDIVIDUALLY PACKAGED

- Single patient use reduces risk of cross-contamination



FOUR INDEPENDENT
SPECIMEN CHAMBERS



DIRECT CONNECTION TO WALL SUCTION
FOR EASY FLUID MANAGEMENT



SEQUENCED TRAP CHAMBERS AND
DIRECT BYPASS CHAMBERS



Non-Sterile, Disposable

Catalog Number	Description	Chambers	Unit
CAP132	Multi Chamber Polyp Trap	4	20/box

Latex **FREE**

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



525 French Road
Utica, NY 13502
(315) 797-8375
Customer Service: 1-800-448-6506
FAX: 1-800-438-3051
International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235
email: info@mail.conmed.com
www.conmed.com