

Rosie™ Scan

MOBILE BLADDER SCANNING TECHNOLOGY

FINALLY! AFFORDABLE, ACCURATE BLADDER SCANNING

This easy to use, non-invasive, mobile touchscreen unit delivers accurate calculation of urine volume to prevent unnecessary catheterizations.

- + Color LED screen displays real-time image
- + Stores 100 cases with patient ID
- + Touchscreen operation
- + Printer on board
- + Accessory tray
- + Rugged rolling cart
- + USB output and storage
- + Rechargeable Lithium-ION battery
- + 2 year warranty on scanner and 12 months on probe

Specifications:

- + Scan Depth: $\geq 140\text{mm}$
- + Sway Angle: $120^\circ \pm 2^\circ$
- + Rotating Angle: $180^\circ \pm 2^\circ$
- + Display Screen: 8.0" Color LCD TFT Screen
- + Patient ID: 12 Digit maximum
- + Case Management: store, browse, support
- + Output: USB
- + Maximum Measurement: 999ml
- + Scan Mode: Male, Female, PRV
- + Transducer Frequency: 2.6 MHz
- + Power Supply: Battery Lithium-ION 4 hours



View a demo of RosieScan on our website or YouTube Channel.

NURSE ROSIE PRODUCTS

DISTRIBUTED BY LIFE SYSTEMS, INC.

7320 CENTRAL AVENUE + SAVANNAH, GA 31406 + 800-841-1109 + NURSE ROSIE.COM