



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Malignant Hyperthermia Cart with Medical Grade Refrigerator #MH5100

Specifications

Construction:

- 20-gauge powder-coated steel with built-in refrigerator

Features:

- Uses high performance Follett® refrigerator or equivalent
- Top drawer allows for storage of poster and MH Procedure Manual available from MHAUS.
- Also store supplies to actively cool a patient, such as hypothermia blankets, over and under the patient
- High performance refrigerator ensures temperature controls for cold isotonic saline for IV infusion and or gastric, peritoneal or rectal irrigation, as appropriate.
- Cart can hold critical MH emergency supplies at the proper temperature in close proximity to areas where they are most needed
- Key lock for drawer
- Cart ships with wheels installed
- Hospital grade plug and cord

Upgrade available:

- Backrail featuring painted steel sign that reads "MH Cart Only"

Drawer configuration:

- (1) – 3" (7.6 cm)
- (1) – Compartment with pre-installed refrigerator

Internal drawer dimensions:

- 23" W x 17" D (58.42 cm W x 43.18 cm D)

Overall cart dimensions:

- 46" H x 30" W x 25"D
- (116.9 cm H x 76.2 cm W x 63.5 cm D)
- Cart net weight: 175 pounds (80 kg)
- Approximate shipping weight: 215 pounds (98 kg)

Please note Refrigerator units are subject to change without notice



Seventeen available color choices:

White, Beige, Sand, Yellow, Light Gray, Mauve,
Light Blue, Teal, Red, Burgundy, Hammer Red, Navy,
Forest Green, HammerTone Gray, HammerTone Green,
HammerTone Blue, Hammer Black