



New DINACLICK Connector!



General-care-reusable BP cuff

GE CRITIKON DURA-CUF
Blood Pressure Cuffs

The first GE CRITIKON* blood pressure cuff, the DURA-CUF* is gentle, yet designed for durability. Made with tough, flexible nylon, they are convenient and easy to clean with curved corners to enhance your patient's comfort. The DURA-CUF is constructed for multiple patient use in low acuity care settings. They are also ideal for the physician office where the risk of contamination is diminished.



Designed for patient safety, comfort and accuracy

- Made with durable materials that can withstand extensive use and multiple disinfections and cleanings
- Designed with rounded corners to enhance patient comfort
- Constructed with materials not made with PVC, making them safer for the environment, and natural rubber latex, making them safer for your patient
- Compliant with AAMI and AHA NIBP cuff sizing standards for prompt, correct sizing

Suitable for virtually every physician's office

- Full range of adapters and connectors are compatible with automated monitors from GE, Philips**, Spacelabs**, and Datascope**, as well as on a variety of manual devices, including wall manometers and hand-held gauges
- Three types of cuff assortment packs enable for easier ordering and optimizing inventory space



DURA-CUF			GE Healthcare				Welch Allyn	Colin Philips/ HP** Spacelabs Siemens** Datascope	Protocol	Manual Sphyg. Wall Manometer
										
Size	Color	Limb Circumference	2-Tube Screw Connector	2-Tube DINACLICK Connector	2-Tube Submin Connector	2-Tube Mated Submin Connector	1-Tube Screw Connector	1-Tube Bayonet Connector	1-Tube Submin Connector	DURA-CUF Inflation System
Infant	Orange	8-13cm	002783	DUR-P1-2A	002200	2750 ¹	002760	002274	2740	002385
Child	Green	12-19 cm	002781	DUR-P2-2A	002201	2751	002761	002275	2741	002384
Child Long	Green	12-19 cm	2785	DUR-P2-2A-L			002769			002386
Small Adult	Royal Blue	17-25 cm	002779	DUR-A1-2A	002202	2752	002762	002276	2742	002383
Small Adult Long	Royal Blue	17-25 cm		DUR-A1-2A-L			002763			002387
Adult w/ Hanger	Navy	23-33 cm	002771							002380
Adult	Navy	23-33 cm	002774	DUR-A2-2A	002203	2753	002764	002277	2743	
Adult Long	Navy	23-33 cm	002772	DUR-A2-2A-L	002206	002756	002765	002130		2388
Large Adult	Wine	31-40 cm	002791	DUR-A3-2A	002204	2754	002766	002278	2744	002381
Large Adult Long	Wine	31-40 cm	002784	DUR-A3-2A-L	002207	002757	002767	002131		2389
Thigh	Brown	38-50 cm	002796	DUR-T1-2A	002205	2755	002768	002279	2745	002382

1. Infant BP Cuff for Datex-Ohmeda monitors has Male Slip Luer Connector

DURA-CUF Assortment Pack	Part Number	Description
2-Tube Screw Connectors	002699	1 Infant, 1 Child, 1 Sm. Adult, 1 Adult, 1 Lg. Adult, and 1 Thigh (6 cuffs/pack)
	002697	2 Infant, 3 Child, 1 Sm. Adult (6 cuffs/pack)
	002698	1 Sm. Adult, 2 Adult, 1 Adult Long, 1 Lg. Adult, 1 Lg. Adult Long (6 cuffs/pack)
2-Tube Submin Connectors	2299	1 Infant, 1 Child, 1 Sm. Adult, 1 Adult, 1 Lg. Adult, and 1 Thigh (6 cuffs/pack)
2-Tube DINACLICK Connectors	2059301-001	1 Sm. Adult, 1 Adult, and 1 Lg. Adult (3 cuffs/pack)

GE Healthcare
P.O. Box 900, FIN-00031 GE, Finland
GE Direct United Kingdom: +44 (0)800 0329201

www.gehealthcare.com



GE imagination at work

To learn more about GE Healthcare's comprehensive cuff solutions, visit us at cuffs.gehealthcare.com

© 2013 General Electric Company — All rights reserved.
General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE representative for the most current information.

GE and GE Monogram are trademarks of General Electric Company.

*Trademarks of General Electric Company.

**All third party trademarks are the property of their respective owner.

EMEA DOC1313654 04/13
(Global version DOC1094825 rev 2)