

Single patient  
short-term use

Single  
patient use

Limited  
reuse

Hospital  
reusable

General care  
reusable



## CRITIKON SENSA-CUF Blood Pressure Cuffs

The CRITIKON\* SENSA-CUF\* is a highly versatile solution for a wide range of extended-care settings. Made using strong, breathable, easy-to-clean material, the SENSA-CUF is designed to maximize patient comfort in low-acuity areas in a hospital where risk of contamination are diminished, and durable enough for multi-patient use.



### A trusted choice for clinical excellence

- Validated to the same quality standard as all GE patient monitors, with repeated inflations to ensure accurate readings
- Large bold print helps you quickly choose the right cuff, enabling a more accurate measurement

### Constructed for patient safety and durability

- Soft, breathable, easy-to-clean material is not made with PVC, making them safer for the environment, or natural rubber latex, making them safer for your patients
- DINACLICK\* connector meets international standards and prevents misconections to luer fittings

### Designed to improve workflow

- Easy-to-use, fast DINACLICK connector enables nurses to link the cuff to the monitor faster than other dual-tube connectors
- Cuffs compliant with AAMI and AHA NIBP cuff sizing standards for prompt, correct sizing
- A full range of adapters and connectors are compatible with automated monitors and a variety of manual devices, including wall manometers, enabling a universal connection system throughout the facility







### Assortment packs available

- For low-volume cuff utilizers, several types of assortment packs are available for easier ordering and more optimal inventory management



DINACLICK  
Connector



SENSA-CUF			GE Healthcare				Non-GE Healthcare	
								
Size	Color	Limb Circumference	2-Tube Screw Connector	2-Tube DINACLICK Connector	2-Tube Submin Connector	2-Tube Mated Submin Connector	1-Tube Screw Connector	1-Tube Bayonet Connector
Infant	Orange	8-13 cm	2458	SEN-P1-2A	2482	2430	2440	2470
Child	Green	12-19 cm	2460	SEN-P2-2A	2484	2431	2441	2472
Small Adult	Lt. Blue	17-25 cm	2462	SEN-A1-2A	2486	2432	2442	2474
Small Adult Long	Lt. Blue	17-25 cm	2463	SEN-A1-2A-L	2487	2433	2443	2475
Adult	Navy	23-33 cm	2464	SEN-A2-2A	2488	2434	2444	2476
Adult Long	Navy	23-33 cm	2465	SEN-A2-2A-L	2489	2435	2445	2477
Large Adult	Rose	31-40 cm	2466	SEN-A3-2A	2490	2436	2446	2478
Large Adult Long	Rose	31-40 cm	2467	SEN-A3-2A-L	2491	2437	2447	2479
Thigh	Brown	38-50 cm	2468	SEN-T1-2A	2492	2438	2448	2480

Note: 5 cuffs/box

SENSA-CUF Assortment Pack	Part Number	Description
2-Tube Screw Connectors	2494	1 Small Adult, 1 Adult, 1 Large Adult (3 cuffs/box)
2-Tube DINACLICK Connectors	2059300-001	1 Small Adult, 1 Adult, 1 Large Adult (3 cuffs/box)

To learn more about GE Healthcare's comprehensive cuff solutions, visit us at [cuffs.gehealthcare.com](http://cuffs.gehealthcare.com)

[www.gehealthcare.com](http://www.gehealthcare.com)

© 2014 General Electric Company — All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your

GE representative for the most current information.

GE and GE Monogram are trademarks of General Electric Company.

\* Trademarks of General Electric Company.

GE Healthcare  
8200 West Tower Avenue  
Milwaukee, WI 53223  
U.S.A.

