

GE SmartMix* for BP Cuffs
Addressing your continuum of patient care



New DINACLICK Connector!



Limited-reuse BP cuff

GE CRITIKON CLASSIC-CUF
Blood Pressure Cuffs

GE CRITIKON* CLASSIC-CUF* BP cuffs are designed and clinically validated for limited reuse on multiple patients. Made with durable vinyl construction, this time-tested CRITIKON blood pressure cuff can also be used on multiple patients. It is both easy to clean and strong to withstand repeated inflations, yet economical enough for single-patient use.



A trusted choice for clinical excellence

- Validated to the same quality standard as all GE monitors with repeated inflations to ensure accurate readings
- Large bold print helps you quickly choose the right cuff, enabling more accurate measurement

Constructed for patient safety and durability


- Easy-to-clean, vinyl material does not contain natural rubber latex or DEHP, making them safer for your patients
- Full range of cuff sizes available, from neonate, pediatric, adult and thigh, ensuring that the measurement you receive is from proper fit
- New DINACLICK* connector meets international standards and prevents misconnections to luer fittings

Designed to improve workflow

- Easy-to-use, fast DINACLICK connector enables nurses to link the cuff to the monitor quicker than a pair of traditional tube connectors
- Cuffs compliant with AAMI and AHA NIBP cuff sizing standards for prompt, correct sizing
- Full range of adapters and connectors are compatible with automated monitors from GE, Philips**, Colin, Protocol, Spacelabs**, Datascope** and Welch Allyn, plus on a variety of manual devices, including wall manometers and hand-held gauges, enables a universal connection system throughout the facility

Cuff assortment packs available

- For low-volume cuff utilizers, several types of assortment packs available for easier ordering and more optimal inventory management



CLASSIC-CUF			GE Healthcare				Welch Allyn	Colin Philips/ HP** Spacelabs Siemens** Datascope	Protocol	Manual Sphyg. Wall Manometer
										
Size	Color	Limb Circumference	2-Tube Screw Connector	2-Tube DINACLICK Connector	2-Tube Submin Connector	2-Tube Mated Submin Connector	1-Tube Screw Connector	1-Tube Bayonet Connector	1-Tube Submin Connector	CLASSIC-CUF Inflation System
Infant	Orange	8-13cm	2618	CLA-P1-2A	2354	2170 ¹	002160	2300	2150	
Child	Green	12-19 cm	2613	CLA-P2-2A	2355	2171	002161	2301	2151	002287
Small Adult	Lt. Blue	17-25 cm	2608	CLA-A1-2A	2356	2172	002162	2302	2152	002285
Adult	Navy	23-33 cm	2603	CLA-A2-2A	2357	2173	002164	2303	2153	002283
Adult Long	Navy	23-33 cm	2647	CLA-A2-2A-L	2352	2174	002165	2306		
Large Adult	Rose	31-40 cm	2643	CLA-A3-2A	2358	2175	002166	2304	2154	002281
Large Adult Long	Rose	31-40 cm	2649	CLA-A3-2A-L	2353	2176	002167	2307		
Thigh	Brown	38-50 cm	2648	CLA-T1-2A	2359	2177	002168	2305	2155	

1. Infant BP Cuff for Datex-Ohmeda monitors has Male Slip Luer Connector

CLASSIC-CUF Assortment Pack	Part Number	Description
2-Tube Screw Connectors	2056304-001	5 Sm Adult, 10 Adult, 5 Lg. Adult (20 cuffs/pack)
	002695	1 Infant, 1 Child, 2 Sm. Adult, 2 Adult, 1 Adult Long, 2 Lg. Adult, 1 Thigh (10 cuff/pack)
2-Tube Submin Connectors	2056304-002	5 Sm Adult, 10 Adult, 5 Lg Adult (20 cuffs/pack)
	2298	1 Infant, 1 Child, 2 Sm. Adult, 2 Adult, 1 Adult Long, 2 Lg. Adult, 1 Thigh (10 cuffs/pack)
2-Tube DINACLICK Connectors	2059304-001	5 Sm. Adult, 10 Adult, 5 Lg. Adult (20 cuffs/pack)

GE Healthcare
P.O. Box 900, FIN-00031 GE, Finland
GE Direct United Kingdom: +44 (0)800 0329201

www.gehealthcare.com



GE imagination at work

To learn more about GE Healthcare's comprehensive cuff solutions, visit us at cuffs.gehealthcare.com

© 2013 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE representative for the most current information.

GE and GE Monogram are trademarks of General Electric Company.

*Trademarks of General Electric Company.

**All third party trademarks are the property of their respective owner.

EMEA DOC1313644 04/13
(Global version DOC1096295 rev 2)