

**ELEVATE**  
The **Health o meter®** EMRscale

#### Includes Connectivity Kit

External adapter allows scale to interface with Welch Allyn Spot Vital Signs® LXi and Welch Allyn Connex® Vital Signs Monitors



US Patent #D667323

### Advantages

#### ✓ EMR Capable

Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions

#### ✓ 350° Swivel Display

Eye-level rotating display can be read from almost any angle and allows for ease of use and patient privacy

#### ✓ Built-in Digital Height Rod

The easy-to-use, state of the art digital height rod provides measurements from 23 5/8" - 88" / 60 cm - 222 cm

#### ✓ Value-Added Features

Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost

#### ✓ Calculates Body Mass Index (BMI)

Scale functions include the ability to calculate the patient's BMI

#### ✓ Patented Design

The ELEVATE-BT's patented design provides a user-friendly and intuitive interface for measuring patient weight, height and BMI with a convenient rotating display head



ELEVATE  
The Health o meter® EMRscale



US Patent #D667323

### Specifications

- Capacity: 600 lb / 272 kg
- Resolution: 0.2 lb / 0.1 kg
- Height Rod Resolution: 0.1" / 0.1 cm
- Platform Size (w x d x h):  
14 1/4" x 14 1/4" x 2 5/8"  
362 mm x 362 mm x 67 mm
- Digital Height Rod: 23 5/8" - 88" / 60 cm - 222 cm
- Includes connectivity kit to interface with Welch Allyn® vital signs monitors
- Connectivity: USB and built-in Bluetooth® wireless technology
- 1" / 25 mm LCD display, 35° left / right swivel  
(350° with height rod fully retracted) for ease of viewing
- Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Hold / Release, Auto Zero, Auto Off
- 2 wheels for ease of movement
- 100 - 240V adapter included, 6 C-cell batteries (not included)
- 2 year limited warranty

Model	Description
SS-ELEVATE-BT	ScaleSurance 2 year extended warranty for ELEVATE-BT models
ADPT50	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
BLUETOOTHONE-P	Wireless thermal printer