



A comprehensive platform scale that supports up to 1000 pounds offers “live” handrails and a low profile platform to provide a stable and comfortable weighing experience for patients of any size. Ideal for any type of facility including primary care, rehab or hospital.

Advantages

- ☒ **EMR Connectivity**
 Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- ☒ **“Live” Handrails**
 Advanced technology allows the patient to hold onto handrails for stability while scale accurately calculates weight
- ☒ **Low Profile Platform**
 Low platform height requires only minimal step up, offering patient a comfortable and stable weighing experience
- ☒ **Calculates Body Mass Index (BMI)**
 Scale functions include the ability to calculate the patient’s BMI
- ☒ **180° Swivel Display**
 The display head has the ability to pivot for the convenience of the caregiver, allowing line-of-sight viewing and information privacy
- ☒ **On-Screen Help Menu**
 Guides users through operating functions and provides easy access to multiple settings
- ☒ **Interactive Display**
 High-contrast color screen provides a user-friendly interface
- ☒ **Value-Added Features**
 Operating functions (LB/KG Conversion, LB/KG LockOut, BMI, Zero,Tare, Pre-tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off, Time/Date, Variable Auto Off Time, Audible/Mute Sound Option) included at no additional cost
- ☒ **Everlock™ (patent pending)**
 Option that allows users to permanently lock the scale in either LB or KG to improve safety



PROPLUSROD:
Optional Height Rod



Specifications

- Capacity: 1000 lb / 454 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size (w x d x h):
22" x 15 3/4" x 2 1/2" / 559 mm x 400 mm x 64 mm
- Product Footprint (w x d x h):
31 1/4" x 27 1/4" x 49 3/4" / 794 mm x 692 mm x 1264 mm
- Product Weight: 52 lb / 23 kg
- Display: 1 1/2" / 38 mm High-Contrast Color TFT-LCD Screen
- Wheels: 2
- Power Source: 100-240V adapter included, 6 D batteries (not included)
- Connectivity: USB, optional Bluetooth®
- Functions: LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off, Time/Date, Variable Auto Off Time, Audible/Mute Sound Option, Everlock™
- 2 year limited warranty

Model	Description
1100KL	Digital platform scale, pounds and kilograms
1100KG	Digital platform scale, kilograms only
1100KLHR	1100KL with optional height rod
1100KGHR	1100KG with optional height rod
SS-1100KL	ScaleSurance 2 year extended warranty (1100 models)*
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
PROPLUSROD	Optional height rod
BLUETOOTHONE-P	Wireless thermal printer