



Advantages

- ✓ **EMR Connectivity**
Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- ✓ **Calculates Body Mass Index (BMI)**
Scale functions include the ability to calculate the patient's Body Mass Index (BMI)
- ✓ **Low Profile Platform**
Low platform height requires only minimal step up, offering patient a more comfortable and stable weighing experience
- ✓ **Value-Added Features**
Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero Out/Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost
- ✓ **Compact Size**
The slim design of the scale and waist-high display require a limited amount of office space



201HR

Optional

Wall Mounted Telescopic
Metal Height Rod



Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
499KL	26.625" x 17.375" x 11.75" / 3.15	22 lb	892076002527



Specifications

- Capacity: 500 lb / 220 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size: 13-3/4" (w) x 16-1/2" (d) x 2-3/8" (h) / 349 mm (w) x 419 mm (d) x 60 mm (h)
- Height Rod: Optional 201HR
- Display: 1" LCD
- Wheels: 2
- Connectivity: USB
- Power Source: 6 AA batteries included, Optional 120V adapter (order # ADPT31)
- Functions: LB / KG Conversion, LB / KG Lock Out, Body Mass Index (BMI), Zero Out / Tare, Hold / Release, Auto Zero, Auto Off
- 2 Year Limited Warranty