



With a high capacity and mechanical height rod, this scale properly weighs and measures the heaviest ambulatory patients. The 350° swivel head makes it easy for clinicians to view measurements and interact with the scale.

Advantages

- ☒ **EMR Capable**
Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions
- ☒ **350° Swivel Display**
Eye-level rotating display can be read from almost any angle and allows for ease of use and patient privacy
- ☒ **Built-in Height Rod**
The easy-to-use height rod provides height measurements from 23 5/8" to 84" / 60cm – 213cm
- ☒ **Value-Added Features**
Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost
- ☒ **Calculates Body Mass Index (BMI)**
Scale functions include the ability to calculate the patient's BMI



350° swivel head display



Specifications



- Capacity: 600 lb / 272 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size (w x h x d):
14 1/4" x 14 1/4" x 2 5/8"
362 mm x 362 mm x 67 mm
- Height Rod: 23 5/8" - 84" / 60 cm - 213 cm
- Display: 1" / 25 mm LCD, 350° Swivel
- Wheels: 2
- Connectivity: USB, optional Bluetooth®
- Power Source: 100-240V adapter included,
6 C-cell batteries (not included)
- Functions: LB / KG Conversion, LB / KG Lock Out , BMI, Zero,
Tare, Hold / Release, Auto Zero, Auto Off
- 2 Year Limited Warranty

Model	Description
SS-597KL	ScaleSurance 2 year extended warranty
ADPT50	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
BLUETOOTHONE-P	Wireless thermal printer