



350° swivel head display



Great for doctor offices and clinics with minimal floor space, the 599KL has a high capacity with waist-high display that can swivel 350°. For providers wishing to measure height, an optional height rod is available.

Advantages

- ☒ **EMR Capable**
 Scale interfaces with a computer, EMR software or other electronic device via Health o meter® Professional Connectivity Solutions
- ☒ **Value-Added Features**
 Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost
- ☒ **350° Swivel Display**
 Eye-level rotating display can be read from almost any angle and allows for ease of use and patient privacy
- ☒ **Calculates Body Mass Index (BMI)**
 Scale functions include the ability to calculate the patient's BMI



Optional height rod for 599KL
(See chart below)



Specifications

- Capacity: 600 lb / 272 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size (w x d x h):
14 1/4" x 14 1/4" x 2 5/8"
362 mm x 362 mm x 67 mm
- Height Rod: Optional 499KLROD
- Display: 1 7/8" / 25 mm LCD, 350° Swivel
- Wheels: 2
- Connectivity: USB, optional Bluetooth®
- Power Source: 100-240V adapter included,
6 C-cell batteries (not included)
- Functions: LB / KG Conversion, LB / KG Lock Out , BMI, Zero, Tare,
Hold / Release, Auto Zero, Auto Off
- 2 year limited warranty

Model	Description
599KLHR	599KL with height rod included
SS-599KL	ScaleSurance 2 year extended warranty
499KLROD	Optional height rod for 599KL
ADPT50	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
BLUETOOTHONE-P	Wireless thermal printer