



A comprehensive platform scale that supports up to 1000 pounds offers “live” handrails and a low profile platform to provide a stable and comfortable weighing experience for patients of any size. Ideal for any type of facility including primary care, rehab or hospital.

Advantages

- ☒ **EMR Connectivity**
 Scale interfaces with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- ☒ **“Live” Handrails**
 Allows the patient to hold onto handrails for stability while scale accurately calculates weight
- ☒ **Swivel Display Head**
 For the convenience of the administrator, the display head rotates and pivots allowing line-of-sight viewing and information privacy
- ☒ **On-Screen Help Menu**
 Guides users through operating functions and provides easy access to multiple settings
- ☒ **Everlock™ (patent pending)**
 Option that allows users to permanently lock the scale in either LB or KG to improve safety
- ☒ **Interactive Display**
 High-contrast color screen provides a user-friendly interface
- ☒ **Value-Added Features**
 Operating functions (LB/KG Conversion, LB/KG Lock Out, (Body Mass Index) BMI, Zero, Tare, Pre-tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off, Time/Date, Variable Auto Off Time, Audible/Mute Sound Option) included at no additional cost



**“LIVE” Handrails
for patient stability**



**Optional height rod
(See chart below)**



Specifications

- Capacity: 1,000 lb / 454 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size (w x d x h):
26" x 22" x 2 3/8" / 660 mm x 559 mm x 60 mm
- Product Footprint (w x d x h):
31" x 27" x 50" / 787 mm x 686 mm x 1270 mm
- Product Weight: 66 lb / 30 kg
- Height Rod: Optional STROD
- Display: 1 1/2" / 25 mm High-Contrast Color TFT-LCD Screen
- Wheels: 2
- Power Source: 100-240V adapter included, 6 D batteries (not included)
- Connectivity: USB, optional Bluetooth®
- Functions: LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off, Time/Date, Variable Auto Off Time, Audible/Mute Sound Option, Everlock™
- 2 Year Limited Warranty

Model	Description
2101KL	Digital platform scale, pounds and kilograms
2101KG	Digital platform scale, kilograms only
2101KLHR	2101KL with height rod included
2101KGHR	2101KG with height rod included
SS-2101KL	ScaleSurance 2 year extended warranty (2101 models)*
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
STROD	Height Rod for 2101KL & 2101KG
BLUETOOTHONE-P	Wireless Thermal Printer