



AirLife® Tri-Flo™ suction catheters and tracheostomy care

AirLife® Tri-Flo™ suction catheters and tracheostomy care kits help you treat the compromised airway. An array of options allows you to choose the optimal product configurations for your patient needs.

AirLife®



Verso™ multi-purpose airway adapter

Help minimize ventilator-associated pneumonia (VAP)

Using the Verso™ adapter in-line allows clinicians to perform open suctioning procedures through a single port *without disconnecting the circuit*.

- Made without DEHP
- Approved for use for up to seven days
- Designed for adult, pediatric and infant use



Adult/Pediatric

Cat. no.	Description	Qty.
CSC100	Verso airway access adapter	20
CSC400	Verso 90 airway access adapter	20
CSC500	Verso T airway access adapter	20

Neonatal/infant

Cat. no.	Description	Qty.
CSC300	Neo-Verso® 'Y' airway access adapter	20
CSC600	Neo-Verso airway access adapter w/swivel	20

Tri-Flo suction catheters

Tri-Flo catheters are designed with two triangular eyes and a beveled tip to increase suctioning efficiency while reducing patient trauma. The Tri-Flo catheter portfolio offers catheters in 5/6, 8, 10, 12, 14 and 18 Fr sizes to cover a range of patient needs. Clear connectors are standard on all 5/6 and 8 Fr sizes to help you easily visualize patient secretions.




Single catheters



Tri-Flo catheters are available in looped and straight packaging configurations in 5/6, 8, 10, 12, 14 and 18 Fr sizes. Most single catheters are available with or without a control port, so you can choose continuous or controlled suctioning.

Looped with control port			Straight with control port		
Cat. no.	Size (Fr)	Qty.	Cat. no.	Size (Fr)	Qty.
T63C	5/6	50	T263C	5/6	100
T64C	8	50	T264C	8	100
T61C	10	50	T261C	10	100
T68C	12	50	T268C	12	100
T60C	14	50	T260C	14	100
T62C	18	50	T262C	18	100
Looped without control port					
T61	10	50			
T60	14	50			
T62	18	50			

No-touch catheters



Tri-Flo catheters are packaged in a plastic sleeve to add a layer of safety between the clinician and the patient. No-touch catheters are available in 8, 10, 12, 14 and 18 Fr sizes and include a pop-up basin for your solution.

With pop-up basin					
Cat. no.	Size (Fr)	Qty.	Cat. no.	Size (Fr)	Qty.
T164C	8	100	T160C	14	100
T161C	10	100	T162C	18	100
T168C	12	100			

Cath-N-Glove™ economy catheter kits



Tri-Flo catheters are packaged with two powder-free vinyl (*latex-free*) exam gloves and a pop-up basin. All economy kits are packaged in a plastic peel pouch to ease use and speed preparation.

Two powder-free vinyl gloves

Cat. no.	Size (Fr)	Qty.
4893T	5/6	100
4897T	8	100
4895T	10	100
4898T	12	100
4894T	14	100

Cath-N-Glove wallet catheter kits



Tri-Flo catheters are packaged with two nitrile (*latex-free*) gloves and a pop-up basin wrapped together in an inner wallet. The kit is packaged in a plastic peel pouch to ease use and speed preparation.

Two nitrile gloves

Cat. no.	Size (Fr)	Qty.
4863T	5/6	100
4867T	8	100
4865T	10	100
4868T	12	100
4864T	14	100

Rigid basin kits



Tri-Flo catheters are packaged in a rigid basin with two walleded powder-free vinyl (*latex-free*) exam gloves. The rigid basin can hold your solution while the inner wallet can provide a sterile field during setup to support aseptic technique.

Two powder-free vinyl gloves

Cat. no.	Size (Fr)	Qty.	Cat. no.	Size (Fr)	Qty.
44-56	5/6	100	44-12	12	100
44-08	8	100	44-14	14	100
44-10	10	100	44-18	18	100

Rigid basin kits with saline



Tri-Flo catheters are packaged in a rigid basin with two walleded powder-free nitrile (*latex-free*) exam gloves and a 4 oz. container of sterile saline. The inner wallet can be used as a sterile field during setup to support aseptic technique.

Two powder-free nitrile gloves

Cat. no.	Size (Fr)	Qty.	Cat. no.	Size (Fr)	Qty.
41-56	5/6	12	41-12	12	12
41-08	8	12	41-14	14	12
41-10	10	12			

Tracheostomy care kits

Complete tracheostomy cleaning trays



We offer a line of sterile, single-use trays, with or without gloves. These trays include:

- Three rigid basins for cleaning, soaking and rinsing
- Pre-sewn trach dressing
- Three gauze sponges
- Three pipe cleaners
- Cleaning brush
- Three cotton-tipped applicators
- Twill tape, 32"
- Plastic forceps
- Water-resistant drape
- Hospital wrap

Cat. no.	Additional components	Qty.
4681A	2 vinyl (<i>latex-free</i>) gloves	20
4682A	Without gloves	20

Water solutions



4 oz foil-lid cups

Cat. no.	Description	Qty.
4872B	Sterile water	48
4873A	Normal saline, 0.9%	48

Tracheostomy care kits




Our kits are supplied in a tray with a coated paper lid. All kits include:

- Trach gauze dressing
- Twill tape, 34"
- Trach brush
- Four gauze sponges, 4" x 4"

Cat. no.	Additional components		Qty.
3T3030A	2 cotton-tipped applicators	3 pipe cleaners	30
3T4691	2 cotton-tipped applicators 2 pipe cleaners	2 vinyl (<i>latex-free</i>) gloves 1 water-resistant towel	20
3T4691A	2 cotton-tipped applicators 3 pipe cleaners	2 vinyl (<i>latex-free</i>) gloves 1 water-resistant towel	30
3T4692	3 cotton-tipped applicators 3 pipe cleaners 2 vinyl (<i>latex-free</i>) gloves	1 water-resistant towel 2 packets hydrogen peroxide, 30 mL 1 trach pad	30

Visit carefusion.com/AirLifeRespiratory to learn more about our unit dose and sterile water and saline solutions.

 **Caution**—U.S. Federal Law restricts these devices to sale by or on the order of a physician.

CareFusion
Vernon Hills, IL

carefusion.com

