

Igloo® Wound Irrigation System



Bionix Igloo Wound Irrigation helps healthcare providers protect staff by preventing back-splash of blood pathogens during high volume irrigation of non-complex wounds.

► PROBLEM:

Simple wounds can become complex problems if harmful bacteria and necrotic tissue are not removed. Wound Irrigation can also put health care workers at risk for cross-contamination of blood-borne pathogens from back-splash.

► SOLUTION:

Igloo Wound Irrigation System utilizes the Igloo Irrigation Shield on a high-volume trigger, which conveniently attaches a variety of saline bottles, allowing you to irrigate at 4-15 psi without having to stop and refill every 60cc's.

► RESULTS:

Continuous high volume wound irrigation that effectively removes bacteria and necrotic tissues, while preventing back-splash and promoting healing.

PRODUCT OFFERING

#5100 – 10 Igloo Wound Irrigation Systems (sterile)

#5250 – 25 Igloo Wound Irrigation Shields (sterile)

#5500 – 50 Igloo Wound Irrigation Shields (sterile)

POTENTIAL REIMBURSEMENT (U.S. ONLY)

12031 – 12057 for intermediate wound repair.

13100 – 13160 for complex wound repair.



*Does Not Include Saline Bottle.

©2016 Bionix Development Corp.



Phone: 800.551.6810
www.Bionix.com

B BIONIX® GOING-BEYOND ►