

# Crosstex® Le Petit Masque

## ***Now There's A Smaller Pink Mask For More Personalized Protection... And A Greater Purpose!***

**The pink Le Petit Masque fits smaller faces perfectly and can even do one more terrific thing — help fight women's cancers!**

The Le Petit Masque is three-quarters of an inch shorter than standard masks to help eliminate gapping and provide a more comfortable, protective fit for smaller faces. And, when you select this Priority Pink product, Crosstex will donate a full 5% of the proceeds to Memorial Sloan-Kettering Cancer Center in support of breast and reproductive cancer research.

### **Le Petit Masque**

- Meets ASTM\* performance classifications as a Low (Primary) barrier medical face mask.
- Meets EN14683 Type II Standard.
- PFE = 98.5% @ 0.1 micron.
- Measures 6.25" x 3.5"/15.875 x 8.89cm compared to standard 7"x3.5"/17.78 x 8.89cm masks.
- Fluid resistant pink outer layer.
- Soft white spunbond inner layer.

### **Ordering Information**

Item #	Quantity	Color
GCAPTPK	50/Box	Pink

\*American Society for Testing and Materials



# Priority Pink

## More Reasons To Purchase Priority Pink... And Continue The Fight Against Women's Cancers!

Your overwhelming support of Priority Pink is helping Crosstex continue to lead the fight against women's cancers. So the more Pink you purchase, the more we can make a difference together in the fight for the cure!

### Pink with a Purpose™ Econoback® Towels

When it comes to Crosstex quality, no details are considered too small.

- Horizontal embossing that allows moisture to flow horizontally across the towel, protecting your patient's clothing.
- Reinforced edge embossing for added strength.
- A unique bonding technique that helps eliminate separation.

Ordering Information		
2 Ply Tissue with Poly: 19 x 13in. 48.26x33.02cm		
Item #	Quantity	Color
WEXPP	500/Case	Pink Ribbons

Choose the Pink Ribbon towel  
and you can now help to wipe  
out women's cancers!



Ordering Information		
Item #	Quantity	Color
ZPP	100/Bag	Pink

### Premium® Pink Saliva Ejectors

Crosstex saliva ejectors are a favorite among dental clinicians, hygienists, and assistants.

- A rounded tip design provides optimal suction without aspirating tissue.
- Hypo-allergenic latex-free construction.
- Soft, non-removable tips are tightly bonded for maximum patient safety.
- Pliable and durable wire-reinforced tubing that is easy to shape, yet maintains the desired form.

Finally, a saliva ejector that's kinder to patients  
and helps you support an important cause!

### Pink with a Purpose™ Sterilization Pouches

Crosstex sterilization pouches set the industry standard for safety, reliability, and compatibility with all steam and EtO sterilization agents.

- Secure seal with triple seal/triple heat processing.
- Wide self-seal adhesive strip.
- Perforated fold.
- Thumb notch.

Ordering Information		
3.5 x 9in./8.89x22.86cm		
Item #	Quantity	Color
SCSPP	200/Box	Pink
5.25 x 10in./13.35x25.4cm		
Item #	Quantity	Color
SCMPP	200/Box	Pink

Keep your practice  
and your patients  
safer while  
supporting a  
worthy cause,  
when you  
purchase Pink!



For more information on the power of ordering Priority Pink,  
visit [www.crosstex.com](http://www.crosstex.com).

10 Ranick Road  
Hauppauge, NY 11788-4209 USA  
Tel: 631-582-6777 • 888-276-7783  
Fax: 631-582-1726 • E-mail: [crosstex@crosstex.com](mailto:crosstex@crosstex.com)

