

Activon Manuka Honey[®]

Wound Healing As Nature Intended



The natural alternative to silver

Distributed by



PASSION. PARTNERSHIP. POSSIBILITIES.

www.dukal.com




www.advancismedical.com

Activon Manuka Honey®

At our state-of-the-art manufacturing facility in the UK, the Advancis range of medical grade Manuka Honey products, offer the unique combination of a high quality performance product at a competitive price.

- ✓ **Effective debridement capabilities for faster wound bed preparation**
- ✓ **Maintains a moist wound healing environment, facilitating faster healing of chronic wounds**
- ✓ **Extensive range of Advancis Manuka honey dressings, allowing the right dressing selection for each wound type**
- ✓ **Honey is non-cytotoxic and is a completely natural substance**

Advancis Medical® bring the benefits of 100% Manuka honey into a versatile suite of wound management tools that enable you to choose the right dressing for every phase of wound healing.



Activon Tulle®

Knitted viscose mesh impregnated with 100% Manuka honey

A knitted viscose mesh primary dressing impregnated with 100% Medical Grade Manuka honey. Activon Tulle® creates a moist healing environment and effectively eliminates wound odor.

Activon Tulle® is ideally suited to debriding and de-sloughing shallow wounds or where the exudate levels have started to decrease.

- ✓ **Activon Tulle® is a 2 in 1 sterile dressing; Manuka honey with a primary dressing**
- ✓ **Allows exudate to pass through the dressing on wounds with low to moderate levels of exudate**
- ✓ **Can be cut to size, increasing patient comfort and resulting in better peri-wound skin condition**
- ✓ **Can be placed on the wound either side up, eliminating incorrect dressing application**



Activon Tube®

100% Medical Grade Manuka honey

Containing 100% Medical Grade Manuka honey with no additives. This wound filler will de-bridge and de-slough, eliminate odors and provides a moist wound healing environment.

Activon Tube® is ideal for debriding necrotic tissue. Activon Tube® can be used as a wound filler for use in cavities, just irrigate with saline solution.

- ✓ **Contains 100% pure Manuka honey with no additives**
- ✓ **Maintains a moist wound healing environment, facilitating faster healing of chronic wounds**
- ✓ **Can be safely used on diabetic patients**
- ✓ **For best results on heavily exuding wounds, use in conjunction with a secondary dressing, in particular the **Eclipse®** super absorbent**

Activon® is a natural alternative to the more common wound therapies such as silver, with none of the perceived challenges that silver products face today, such as their cost-effectiveness, lack of efficacy and concerns around their safety and resistance¹.

Item #	Description	Packaging	HCPCS
CR3830	Activon Tube®, 0.9 oz.	12/bx	A9270
CR4135	Activon Tulle®, 2" x 2"	5/bx	A9270
CR4136	Activon Tulle®, 4" x 4"	5/bx	A9270

¹ Cooper Rose; Gray David. 'Is Manuka Honey a Credible Alternative to Silver in Wound Care?' 2012, Wounds UK.

Your patient outcomes can be improved when our Activon Manuka Honey® dressings are used in conjunction with our **Eclipse® range of super-absorbers.**



DUKAL Corporation
2 Fleetwood Court
Ronkonkoma, NY 11779
Tel: 1-800-243-0741 | Email: advancisinfo@dukal.com

MAR436

Advancis
Medical
www.advancismedical.com