

DEVRASORB™

Foam Dressings



Features & Benefits

- Super absorbency – absorbs more than 15 times and retains more than 10 times its weight in fluids
- Optimized wound breathability – outer semipermeable membrane limits liquid leakage yet allows physiologic oxygen and moisture vapor permeability
- Wound cleansing – incorporates a nonionic surfactant that has been found to be an effective wound cleansing agent
- Ease of release – provides painless release even after several days
- Contains skin soothing agents – skin softeners and humectants exclusively formulated to prevent the dressing from drying out and sticking to the wound bed
- Isotonic solutions create a safe, physiologic environment that has long been known to be beneficial to cellular healing

* DevraSorb™ is a registered trademark of Curaline Inc.

Item #	Description	Packaging	HCPCS
Sterile - Latex Free			
1010303	Non-Adhesive Pad, 3" x 3"	10/bx, 5bx/cs	A6209
1010404	Non-Adhesive Pad, 4" x 4"	10/bx, 5bx/cs	A6209
1010405	Non-Adhesive Pad, 4" x 5"	10/bx, 5bx/cs	A6210
1010413	Non-Adhesive Pad, 4" x 13"	10/bx, 5bx/cs	A6211
1010606	Non-Adhesive Pad, 6" x 6"	10/bx, 5bx/cs	A6210
1022525	Adhesive Island, 4" x 4"	10/bx, 5bx/cs	A6212
1024545	Adhesive Island, 6" x 6"	10/bx, 5bx/cs	A6213
2010405	Ultra Non-Adhesive Pad, 4" x 5"	10/bx, 5bx/cs	A6210
2010606	Ultra Non-Adhesive Pad, 6" x 6"	10/bx, 5bx/cs	A6210
3010108	Prime Non-Adhesive Pad, 1" x 8"	10/bx, 5bx/cs	A6215

Connect with us!

Find our full line of products online at:
www.dukal.com

Telephone: 1-631-656-3800

Toll Free: 1-800-243-0741

Fax: 1-631-656-3810

Email: customerservice@dukal.com

To Place Orders: orders@dukal.com

Corporate Office

DUKAL Corporation
2 Fleetwood Court,
Ronkonkoma, NY 11779



STRATEGIC
Product Partnership