



# FLEXAFORM™-ES ELASTIC TAPE



Flexaform™-ES Tape is made of a durable and porous cotton elastic with non-latex adhesive. Ideal for a variety of athletic and orthopedic uses, Flexaform™-ES is highly conformable with superior strength and can be easily removed.

## **Non-Sterile - Latex Free**

ITEM #	DESCRIPTION	PACKAGING
8211	1" x 2.5 yds (5yds stretched), tan	12 rl/bx, 12 bx/cs
8212	2" x 2.5 yds (5yds stretched), tan	6 rl/bx, 12 bx/cs
8213	3" x 2.5 yds (5yds stretched), tan	4 rl/bx, 12 bx/cs
8214	4" x 2.5 yds (5yds stretched), tan	6 rl/bx, 6 bx/cs

## **Features & Benefits:**

- Superior strength & controlled compression
- High conformability for all body contours
- Full unwind allows for continually consistent taping
- Porous construction permits moisture transmission & sustains skin integrity
- Red line centered on the back of cloth provides guide for uniform application and pressure
- Adhesive is strong but doesn't restrict movement
- Non-Latex adhesive reduces skin irritation

## **Suggested Markets & Uses:**

- Orthopedic related injuries
- Physical Therapy & Sports Medicine
- Added athletic support
- Optimal support for fractures, muscle injuries, knees, elbows, and shoulders
- Veterinary Care

Quality



Value



Performance

