



PASSION. PARTNERSHIP. POSSIBILITIES.



# Wavelength<sup>®</sup> Ultrasound Gel

## Description

Wavelength<sup>®</sup> Ultrasound Gel is an excellent conductive agent for use in therapeutic ultrasound, diagnostic imaging, and other external procedures that require a thick water soluble medium.

## Features & Benefits

- Salt free, Fragrance free, Alcohol free, Formaldehyde free
- Water soluble, hypoallergenic, bacteriostatic
- Conductive, non-corrosive, greaseless, and non-staining
- Safe and effective for all transducers
- Flip top cap with silicone barrier to prevent drawback after pressure is released from bottle
- Thick, viscous gel that spreads evenly without running or drying
- Available in blue or clear
- Approved by major diagnostic ultrasound equipment manufacturers

## Contraindications

- Not to be used for an electromedical conductor
- Not to be used for invasive ultrasound procedures

Item	Description	Packaging
Non-Sterile - Latex Free		
NTPA201W	Blue, Flip Cap, Linerless 250ml, 8.5 oz.	12/cs
NTPC201W	Clear, Bottle, Flip Cap, Linerless, 250ml, 8.5 oz.	12/cs
NTPA025X	Blue Single Use Packet 25 g	100/cs
Non-Sterile - Latex Free - Master Packed; w/4 Empty Dispenser Bottles		
NTPA503X	Blue, Cube, 5L	4/cs
NTPC503X	Clear, Cube, 5L	4/cs
NTPF503X	Clear, High Viscosity Gel, Cube, 5L	4/cs
Non-Sterile - Latex Free - Master Packed; w/4 Empty Dispenser Bottles		
NTPD503X	Ultrasound Lotion, 5L	4/cs
Non-Sterile - Latex Free		
NTPD352X	Ultrasound Lotion, 1GAL	4/cs

## Connect with us!

Find our full line of products online at:  
[www.dukal.com](http://www.dukal.com)

Telephone: 1-631-656-3800

Toll Free: 1-800-243-0741

Fax: 1-631-656-3810

Email: [customerservice@dukal.com](mailto:customerservice@dukal.com)

To Place Orders: [orders@dukal.com](mailto:orders@dukal.com)

## Corporate Office

DUKAL Corporation  
2 Fleetwood Court  
Ronkonkoma, NY 11779



**STRATEGIC**  
Product Partnership